

Memmert U Universal Ovens Continued



Accessories

Further options and accessories are available on request. Additional equipment must be specified when ordering.

Code	Alt Ref	Description	For use with	Price
INC8202	E28884	Stainless steel grid, electropolished	30 models	£95.00
CHA4020	E20164	Stainless steel grid, electropolished	55, 75 models	£100.00
CHA4022	E20165	Stainless steel grid, electropolished	110, 160 models	£134.00
CHA4052	E28891	Stainless steel grid, electropolished	260 models	£160.00
CHA4054	E20182	Stainless steel grid, electropolished	450, 750 models	£246.00
INC8128	B03916	Perforated stainless steel shelf	55, 75 models	£103.00
CHA4062	B29725	Perforated stainless steel shelf	260 models	£171.00
CHA4064	B00328	Perforated stainless steel shelf	450, 750 models	£261.00

Memmert VO Vacuum Ovens



Reliable, comfortable and energy efficient - that's how vacuum drying works today. Separately adjustable and directly heated thermoshelves ensure precise temperature control and, thanks to the digital pressure control, the chamber load is not only dried at high speed but also extremely gently in the powerful and robust vacuum drying oven. Together with the speed-controlled vacuum pump (sold separately), the vacuum oven VO forms an unbeatably energy-efficient duo.

- Temperature range: up to 200°C
- Vacuum control range: 5 to 1100 mbar
- 1 model variant: TwinDISPLAY
- Anti splinter; VDE-tested door construction for all models
- Pump control: optimised rinsing of the pump membrane as well as signal output for switching the pump ON/OFF according to requirements
- User friendly, intuitive operating concept
- Nearly exclusive use of high quality, rust resistant and easy to clean stainless steel for interior and exterior housing
- Precise and homogenous temperature control thanks to a product-specific heating concept
- Wide range of options for programming and documentation using interfaces, integrated data loggers and the AtmoCONTROL software
- Cooled vacuum oven VOcool also available (please contact customer services for further information)
- Further options and accessories are available on request, additional equipment must be specified when ordering



For further information on TwinDISPLAY please see page 529

Code	Alt Ref	Volume, L	Int Dims, w x d x h, mm	Ext Dims, w x d x h, mm	Price
OVE1100	VO29	29	385 x 250 x 305	550 x 400 x 607	£6660.00
OVE1102	VO49	49	385 x 330 x 385	550 x 480 x 687	£7520.00
OVE1104	VO101	101	545 x 400 x 465	710 x 550 x 767	£9992.00