

## Genlab Large Capacity Ovens



The Genlab large capacity ovens offers a selection of highly efficient, reliable and cost effective ovens to suit most warming, drying and heat treatment processes. The control system comprises of a touch screen interface that offers both accurate temperature control and integral overheat system by using 2 individual Pt100 sensors. Advance features can be turned on and off to suit bespoke heating applications and tailored to customer preference. The vertical style units are fitted with single doors, the horizontal style units are fitted with twin doors.

- Temperature range: 50 to 250°C
- On screen historical trending (48 hours)
- 8 stage profile control
- Countdown timer with 24/7 real time clock operation
- Manual or automatic overheat reset
- High accuracy Pt100B duplex sensors <0.8°C
- Fan assisted circulation
- Excellent stability <±0.6°C
- Low chamber uniformity
- Energy efficient with integral energy meter
- New - Set point lock feature



**Vertical Style with Single Door**

Code	Alt Ref	Capacity, L	No. of Shelves	Shelf Positions	Int Dims, w x d x h, mm	Ext Dims, w x d x h, mm	Price
<b>OVE2026</b>	LCQ/250V/TDIG	240	5	7	510 x 520 x 910	980 x 690 x 1170	<b>£3313.00</b>
<b>OVE2028</b>	LCQ/350V/TDIG	345	5	7	610 x 620 x 910	1080 x 790 x 1170	<b>£3575.00</b>
<b>OVE2030</b>	LCQ/425V/TDIG	425	5	7	610 x 770 x 910	1130 x 950 x 1270	<b>£4542.00</b>
<b>OVE2032</b>	LCQ/500V/TDIG	535	5	9	610 x 770 x 1140	1130 x 950 x 1400	<b>£4891.00</b>
<b>OVE2034</b>	LCQ/650V/TDIG	675	5	9	770 x 770 x 1140	1450 x 950 x 1400	<b>£6869.00</b>
<b>OVE2036</b>	LCQ/750V/TDIG	760	5	10	770 x 770 x 1280	1450 x 950 x 1640	<b>£7091.00</b>
<b>OVE2038</b>	LCQ/850V/TDIG	840	5	11	770 x 770 x 1420	1450 x 950 x 1780	<b>£7630.00</b>
<b>OVE2040</b>	LCQ/1000V/TDIG	1005	5	11	770 x 920 x 1420	1450 x 1100 x 1780	<b>£8160.00</b>

## Chemicals & Consumables



Laboratory consumables  
and chemicals to  
suit a wide range  
of applications