

Genlab Dual Purpose Ovens/Incubators Continued



Side Mount Design

198L version includes fan assisted circulation as standard.



Code	Alt Ref	Capacity, L	No. of Shelves	Shelf Positions	Ext Dims, w x d x h, mm	Int Dims, w x d x h, mm	Price
INC2280	DP/50/TDIG	58	2	3	730 x 480 x 580	400 x 350 x 420	£1290.00
INC2282	DP/75/TDIG	78	2	3	850 x 480 x 580	530 x 350 x 420	£1496.00
INC2284	DP/100/TDIG	102	3	4	850 x 590 x 580	530 x 460 x 420	£1635.00
INC2286	DP/125/TDIG	126	3	5	850 x 590 x 690	530 x 460 x 520	£1773.00
INC2288	DP/150/TDIG	149	3	5	850 x 670 x 690	530 x 540 x 520	£1968.00
INC2290	DP/200/F/TDIG	198	4	7	850 x 670 x 890	490 x 540 x 750	£2485.00



Genlab Prime Ovens



The flagship brand of Genlab ovens, built to the highest standards and packed with many features and enhancements. As standard, the Prime range offers a heavy duty stainless steel interior, increased accuracies on the standard models and are quicker responding. This range can offer increased load capability, multiple reinforced shelves and operate at 250/300°C continuously. It features our latest touch screen control system which offers intuitive control, designed for industrial and laboratory ovens.

- Stainless steel internal chamber
- Temperature range: choice of 40 to 250°C or 40 to 300°C
- On screen historical trending (48 hours)
- Countdown timer
- 24/7 real time clock operation
- Manual or automatic overheat reset
- High accuracy Pt100A duplex sensor <0.6°C
- Fan assisted circulation
- Excellent stability <±0.5°C
- Low chamber uniformity
- Quick responding
- Energy efficient with integral energy meter
- Set point lock feature

Code	Alt Ref	Temp Range, °C	Capacity, L	No. of Shelves	Shelf Positions	Ext Dims, w x d x h, mm	Int Dims, w x d x h, mm	Price
OVE2150	PRO/30/TDIG	40 - 250	30	2	2	500 x 620 x 570	360 x 350 x 240	£1374.00
OVE2152	PRO/50/TDIG	40 - 250	53	2	3	630 x 590 x 660	490 x 330 x 330	£1552.00
OVE2154	PRO/75/TDIG	40 - 250	73	2	3	630 x 720 x 660	490 x 450 x 330	£1707.00
OVE2156	PRO/100/TDIG	40 - 250	99	3	4	630 x 720 x 780	490 x 450 x 450	£1807.00
OVE2158	PRO/125/TDIG	40 - 250	120	3	5	630 x 720 x 880	490 x 450 x 550	£1963.00
OVE2160	PRO/150/TDIG	40 - 250	143	3	5	630 x 800 x 880	490 x 530 x 550	£2124.00
OVE2162	PRO/200/TDIG	40 - 250	195	4	7	630 x 800 x 1080	490 x 530 x 750	£2529.00
OVE2084	PRO/30/TDIG/300	40 - 300	30	2	2	500 x 620 x 570	360 x 350 x 240	£1769.00
OVE2086	PRO/50/TDIG/300	40 - 300	53	2	3	630 x 590 x 660	490 x 330 x 330	£1945.00
OVE2088	PRO/75/TDIG/300	40 - 300	73	2	3	630 x 720 x 660	490 x 450 x 330	£2101.00
OVE2090	PRO/100/TDIG/300	40 - 300	99	3	4	630 x 720 x 780	490 x 450 x 450	£2202.00
OVE2092	PRO/125/TDIG/300	40 - 300	120	3	5	630 x 720 x 880	490 x 450 x 550	£2356.00
OVE2094	PRO/150/TDIG/300	40 - 300	143	3	5	630 x 800 x 880	490 x 530 x 550	£2519.00
OVE2096	PRO/200/TDIG/300	40 - 300	195	4	7	630 x 800 x 1080	490 x 530 x 750	£3058.00