


Meridian Bioscience ISOLATE II Plant DNA Kits

meridian BIOSCIENCE[®]
A LIFECORE TECHNOLOGIES COMPANY

ISOLATE II Plant DNA Kit provides a simple, efficient column-based method for the isolation of genomic DNA from a wide variety of plant materials, without the need for hazardous reagents such as phenol

- Fast - streamlined protocol for reliable recovery of genomic DNA in as little as 30 minutes
- High-performance - extraction of high-quality DNA, ideal for use in all downstream applications
- Efficient - optimised lysis conditions and column matrix for improved recovery of genomic DNA from a wide range of plant samples
- Versatile - a choice of two lysis buffers ensures excellent genomic DNA recovery from any plant material
- Convenient - kit includes all necessary components, including filters (shredders) and RNase A
- Safe - no hazardous phenol/chloroform extraction




Code	Pack	Price
BI052068	10 preps	£47.00
BI052069	50 preps	£245.00
BI052070	250 preps	£1113.00

Meridian Bioscience ISOLATE II RNA Mini Kits

meridian BIOSCIENCE[®]
A LIFECORE TECHNOLOGIES COMPANY

Highly efficient purification of total RNA by combining the stringency of guanidinium thiocyanate lysis with the speed and purity of silica-membrane column purification.

- Fast - simple extraction of high purity total RNA from six samples, ideal for use in all applications, in as little as 30 minutes
- High-performance - recovery of consistently high quality RNA from a wide variety of animal, bacterial and plant cells and tissues and cell-free samples
- Efficient - optimised lysis conditions and column matrix for recovery of up to 70mg of high purity RNA from each sample
- Convenient - includes all necessary components, including filters (shredders) and DNase I
- Safe - no hazardous phenol/chloroform extraction, CsCl centrifugation, LiCl or alcohol precipitation



Code	Pack	Price
BI052071	10 preps	£69.00
BI052072	50 preps	£320.00
BI052073	250 preps	£1447.00

Electrophoresis

A complete range of electrophoresis tanks, power supplies and accessories brought to you by SLS.

See pages **341-382**

