

Eppendorf twin.tec® PCR Plates

eppendorf

These plates combine very thin walled polypropylene wells for optimal heat transfer to the sample with robust polycarbonate frames for ultimate rigidity and torque-resistance. Twin.tec® plates reduce the risk of cross contamination from their raised well rims for effective sealing.

- One-piece design for optimum performance
- Minimal binding of DNA, RNA and enzymes
- Cut away corner and alphanumeric grid referencing
- Autoclavable (121°C, 20 min)
- Free of detectable human DNA, DNase, RNase and PCR inhibitors
- Skirted plates are compatible with automated systems, thermal cyclers and MegaBACE™ capillary sequence
- Semi-skirted plates are ideal for qRT PCR and fit most thermal cyclers
- Microbiology PCR plates are supplied sterile, individually blistered, PCR clean and free of any bacterial DNA
- 384-well plates have a U-shaped inner contour of the well for enhanced sample visibility



Code	No. of Wells	Description	Well Volume, µL	Colour	Pack	Price
E0030128648	96	Skirted	150	Colourless	25	£119.00
E0030128656	96	Skirted	150	Yellow	25	£119.00
E0030128664	96	Skirted	150	Green	25	£119.00
E0030128672	96	Skirted	150	Blue	25	£119.00
E0030128680	96	Skirted	150	Red	25	£119.00
E0030129300	96	Skirted, microbiology	150	Colourless	10	£72.00
E0030129318	96	Skirted, microbiology	150	Blue	10	£72.00
E0030128575	96	Semi-skirted	250	Colourless	25	£119.00
E0030128583	96	Semi-skirted	250	Yellow	25	£119.00
E0030128591	96	Semi-skirted	250	Green	25	£119.00
E0030128605	96	Semi-skirted	250	Blue	25	£119.00
E0030128613	96	Semi-skirted	250	Red	25	£119.00
E0030129326	96	Semi-skirted, microbiology	250	Colourless	10	£72.00
E0030129334	96	Semi-skirted, microbiology	250	Blue	10	£72.00
E0030133366	96	Unskirted	250	Colourless	20	£76.00
E0030133390	96	Unskirted	250	Blue	20	£76.00
E0030133307	96	Unskirted, low profile	150	Colourless	20	£76.00
E0030133331	96	Unskirted, low profile	150	Blue	20	£76.00
E0030133374	96	Unskirted, divisible	250	Colourless	20	£136.00
E0030133404	96	Unskirted, divisible	250	Blue	20	£136.00
E0030133358	96	Unskirted, low profile, divisible	150	Colourless	20	£136.00
E0030133382	96	Unskirted, low profile, divisible	150	Blue	20	£136.00
E0030128508	384	Skirted	40	Colourless	25	£183.00
E0030128516	384	Skirted	40	Yellow	25	£183.00
E0030128524	384	Skirted	40	Green	25	£183.00
E0030128532	384	Skirted	40	Blue	25	£183.00
E0030128540	384	Skirted	40	Red	25	£183.00
E0030129350	384	Skirted, microbiology	40	Blue	10	£106.00