

Memmert WTB Unstirred Water Baths

A laboratory water bath with the highest degree of comfort. Perfect for everyday use in the laboratory: the high quality stainless steel container with rounded corners is easy to clean and can be completely drained using the standard drain valve without residue left behind. The compact ergonomic design means the Memmert WTB fits on any laboratory table. The touchscreen impresses with its intuitive and convenient menu navigation. The modular concept of this modern water bath guarantees maximum flexibility. If the application calls for it, a circulation pump, cooling unit, shaking device and test tube holder or insertable rack can easily be retrofitted.

- Temperature range up to +100°C
- State of the art control technology
- High quality, corrosion resistant and easily cleanable stainless steel for the working chamber and housing
- Extensive range of accessories sold separately: e.g. a Peltier cooling unit, flat or gable cover, shaking device, circulation pump



Code	Alt Ref	Volume, L	Ext Dims, w x d x h, mm	Int Dims, w x d x h, mm	Price
BAT8300	WTB6	7	349 x 289 x 328	227 x 217 x 150	£896.00
BAT8302	WTB11	10	349 x 289 x 398	267 x 207 x 200	£946.00
BAT8304	WTB15	17	554 x 349 x 348	482 x 277 x 150	£1018.00
BAT8306	WTB24	23	554 x 349 x 398	472 x 267 x 200	£1219.00
BAT8308	WTB35	37	674 x 554 x 348	602 x 482 x 150	£1340.00
BAT8310	WTB50	51	674 x 554 x 398	592 x 472 x 200	£1456.00

SLS Lab Pro Stirred Water Baths

Featuring digital PID temperature control. The thermostirrer is robustly constructed from stainless steel for durability. Painted surfaces feature a tough antibacterial finish with good chemical resistance and immersed components are corrosion resistant. Quiet running in use and providing water temperature stability in the bath, an alarm indicates over temperature and timer running. Low water level is indicated and heating is automatically suspended with an automatic restart when the bath is filled.

- Sensitivity: $\pm 0.05^{\circ}\text{C}$
- Uniformity: $\pm 0.01^{\circ}\text{C}$
- Temperature range: ambient +5 to 99°C
- LED splash proof display with PID temperature control
- Can be set for water temperature or a combination of temperature and time
- Timer can be programmed to commence once a preset temperature is reached and setting retained once unit is turned off
- Under/over temperature alarm is automatically set 4°C above/below set point
- Includes a stainless steel shelf ensuring good water circulation and temperature stability



Code	Capacity, L	Ext Dims, w x d x h, mm	Int Dims, w x d x h, mm	Price
SLS1170	8	271 x 332 x 170	129 x 300 x 150	£891.00
SLS1172	14	361 x 332 x 170	219 x 298 x 150	£914.00
SLS1174	22	537 x 332 x 170	395 x 298 x 150	£963.00

Lids

Lift off, stainless steel lids.

Code	For use with	Price
SLS1186	8L bath	£143.00
SLS1188	14L bath	£157.00
SLS1190	22L bath	£166.00