

Eppendorf HeatSealer Microplate Sealer

eppendorf

The Eppendorf HeatSealer Family offers you reliable and effective sealing of a variety of plates. You can choose from a more compact (S100) and a more flexible (S200) instrument according to your sealing needs.

- Hermetic sealing of multi-well plates
- Minimises evaporation in PCR, reducing cross contamination
- Seals 96- and 384-well plates
- Integrated thermostat prevents overheating
- Improved mechanics for easier sealing
- PCR plate adapter included



Code	Model	Temperature, °C	Weight, kg	Dims, w x d x h, mm	Power, W	Price
E5391000036	S100	170	4	220 x 250 x 230	320	£2302.00
E5392000030	S200	125 to 200	7.2	220 x 321 x 425	350	£3304.00

Accessories

Code	Description	For use with	Price
E5391070018	Low profile plate adapter	S100	£294.00
E5392070020	Low profile plate adapter	S200	£294.00
E5391070034	High profile plate adapter	S100	£294.00
E5392070011	High profile plate adapter	S200	£294.00
E5391070026	96-well PCR plate adapter	S100	£294.00
E5392070038	96-well PCR plate adapter	S200	£294.00


Eppendorf Heat Sealing Foil/Film

eppendorf

- Hermetic sealing of multi-well plates, especially recommended for low reaction volumes
- Best protection against evaporation during PCR
- Maximum adhesiveness at high temperatures during PCR

Supplied in packs of 100

Code	Description	Price
E0030127854	Heat sealing foil	£114.00
E0030127838	Heat sealing film	£162.00



BRAND 0.2mL PCR Tubes



- Free of DNase, DNA, RNase, and endotoxins for reliable and contaminant free results
- Thin walls for rapid heat transfer
- Tightly closing covers protect samples against evaporation

Supplied in packs of 1000

Code	Alt Ref	Description	Price
MOL2300	781300	Domed cap	£73.00
MOL2302	781305	Flat cap	£60.00

