




Azenta Life Sciences Strips of Sealing Caps



- Strips of 8 or 12 domed or flat optical caps
- Universal fit for all 4titude® PCR tube strips and PCR plates
- Crystal clear strips are made of a special polymer with improved optical properties and evaporation-safe fit leading to high transmission rates; ideally suited for small samples with low signal intensity


Code	Alt Ref	Description	Pack	Price
PCR1066	4ti-0751	Strips of 8 flat optical caps	300	£89.00
PCR1176	4ti-0788	Strips of 12 flat optical caps	200	£60.00
PCR1068	4ti-0752	Strips of 8 domed caps	300	£89.00
PCR1092	4ti-0755	Crystal clear strips of 8 flat optical caps	300	£93.00
PCR1094	4ti-0755/120	Crystal clear strips of 8 flat optical caps	120	£41.40



Azenta Life Sciences FrameStar® 96 Well Semi-Skirted PCR Plates, ABI® Style



- Specifically designed to be directly compatible with all major thermal cyclers, this plate can be used directly in ABI® 96 well instruments without the need for any adapters
- Rigid FrameStar® 2 component design eliminates warping and distortion during PCR
- Semi-skirt allows for labeling or barcoding for sample tracking
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen
- Ultra smooth, uniform, thin walled polypropylene
- Standard profile, cut corner at A12 with rigid polycarbonate frame
- Reduced thermal expansion and sample evaporation for improved consistency
- Improved seal integrity due to thermal stability of frame
- Alphanumeric grid reference to aid well and sample identification
- Working volume: <200µL
- Total well capacity: 300µL



Code	Alt Ref	Well Colour	Frame Colour	Pack	Price
PCR1104	4ti-0770/C/10P	Clear	Clear	10	£30.90
PCR1102	4ti-0770/C	Clear	Clear	50	£219.00
PCR1118	4ti-0771	White	Black	50	£150.00
PCR1122	4ti-0772/10P	Frosted	Clear	10	£69.00
PCR1120	4ti-0772	Frosted	Clear	50	£220.00



Azenta Life Sciences 96 Well Semi-Skirted PCR Plates



- Working volume: <200µL
- Total well capacity: 300µL
- Universal semi-skirted plate
- Designed for use on standard thermal cyclers
- Manufactured from virgin polypropylene in Class 7 ISO certified clean room
- Certified RNase, DNase and human genomic DNA free

Supplied in packs of 50

Code	Alt Ref	Colour	Price
PCR1096	4ti-0760	Clear	£99.00
PCR1098	4ti-0761	White	£111.00