


Azenta Life Sciences 96 Well Non-Skirted PCR Plates



- 96 well non-skirted PCR plate, standard profile, 0.2mL wells, polypropylene, cut corner H12
- Working volume: <200µL
- Total well capacity: 300µL
- Universal non-skirted plate
- Designed for use on standard thermal cyclers
- Manufactured from virgin polypropylene in Class 7 ISO certified clean room
- Certified RNase, DNase and human genomic DNA free

Code	Alt Ref	Colour	Pack	Price
PCR0984	4Ti-0750	Clear	50	£99.00



Azenta Life Sciences FrameStar® 96 Well Skirted PCR Plates



- Low profile wells are shorter than "standard" profile wells, which decrease the "dead space" between the heated lid of the thermal cycler and the sample in the well, this eliminates condensation forming on the side wall of the tubes, preventing reduction in PCR volume and increasing the efficiency of the reaction
- Especially recommended for reaction volumes below 20µL
- Rigid polycarbonate skirt of this design eliminates warping and distortion during PCR
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen
- Stackable
- Ultra-smooth, uniform, thin-walled polypropylene
- Working volume: <100µL
- Total well capacity: 200µL
- Rigid polycarbonate frame for added mechanical stability, cut corner H1
- Significantly reduced thermal expansion and sample evaporation for improved consistency in PCR results
- Improved seal integrity due to thermal stability of frame
- Alphanumeric grid reference to aid well and sample identification



Supplied in packs of 50

Code	Alt Ref	Description	Well Colour	Frame Colour	Price
PCR1220	4ti-0960	Fully skirted	Clear	Purple	£145.00
PCR1232	4ti-0960/C	Fully skirted	Clear	Clear	£219.00
PCR1234	4ti-0960/RIG	Fully skirted with rigid skirt	Clear	Black	£229.00
PCR1236	4ti-LB0960/RIG	Fully skirted extra rigid low binding	Clear	Black	£263.00
PCR1244	4ti-0961	Fully skirted	White	Black	£246.00
PCR1246	4ti-0966	Fully skirted	Black	Black	£147.00

Azenta Life Sciences 96 Well Skirted PCR Plates



- Working volume: <100µL
- Total well capacity: 200µL
- Universal 96 well skirted plate
- Designed for use on standard thermal cyclers
- Manufactured from virgin polypropylene in Class 7 ISO certified clean room
- Certified RNase, DNase and human genomic DNA free

Supplied in packs of 50

Code	Alt Ref	Colour	Price
PCR0978	4ti-0740	Clear	£99.00
PCR0980	4ti-0741	White	£111.00

