


## MBL® Class B Pipettes



- Complies with ISO 835
- Graduations in permanent amber stain
- Manufactured from soda-lime glass with robust tooled jets



### Type 1

- Calibrated to deliver from zero at the top to any graduation line down to the shoulder

#### Supplied in packs of 5

Code	Alt Ref	Capacity, mL	Graduations, mL	Tolerance, ± mL	Colour Code	Price
PIP22000	PRA760B	1	0.01	0.01	Yellow	£13.60
PIP22020	PRA762B	2	0.02	0.02	Black	£15.40
PIP22040	PRA764B	5	0.05	0.05	Red	£15.80
PIP22060	PRA766B	10	0.10	0.10	Orange	£16.60
PIP22080	PRA768B	25	0.10	0.20	Green	£24.40

### Type 2

- Calibrated to deliver from any graduation line down to zero at the jet

#### Supplied in packs of 5

Code	Alt Ref	Capacity, mL	Graduations, mL	Tolerance, ± mL	Colour Code	Price
PIP23100	PRA800B	1	0.01	0.01	Yellow	£13.60
PIP23120	PRA802B	2	0.02	0.02	Black	£15.40
PIP23140	PRA804B	5	0.05	0.05	Red	£15.80
PIP23160	PRA806B	10	0.10	0.10	Orange	£16.60
PIP23180	PRA808B	25	0.10	0.20	Green	£24.40

### Type 3

- Calibrated to deliver from zero at the top to any graduation line down to the jet

#### Supplied in packs of 5

Code	Alt Ref	Capacity, mL	Graduations, mL	Tolerance, ± mL	Colour Code	Price
PIP23200	PRA830B	1	0.01	0.01	Yellow	£16.80
PIP23220	PRA832B	2	0.02	0.02	Black	£16.80
PIP23240	PRA834B	5	0.05	0.05	Red	£16.80
PIP23260	PRA836B	10	0.10	0.10	Orange	£18.50
PIP23280	PRA838B	25	0.10	0.20	Green	£30.30

### Type 4, Serological

- Calibrated to deliver from any graduation line down to zero at the jet with the last drop expelled by blowing
- Tooled suction ends accept cotton-wool plug

#### Supplied in packs of 5

Code	Alt Ref	Capacity, mL	Graduations, mL	Tolerance, ± mL	Colour Code	Price
PIP23400	PRA840S	1	0.01	0.01	Yellow	£16.60
PIP23420	PRA842S	2	0.02	0.02	Black	£17.80
PIP23440	PRA844S	5	0.05	0.05	Red	£19.70
PIP23460	PRA846S	10	0.10	0.10	Orange	£20.40
PIP23480	PRA848S	25	0.10	0.20	Green	£36.40