

Grant Optima Stirred Water Baths Continued



Optima Tanks

Baths with intelligent design, robust construction and an excellent working area. Available in stainless steel or plastic.



Code	Alt Ref	Capacity, L	Material	Ext Dims, w x d x h, mm	Int Dims, w x d x h, mm	Drain Outlet	Price
BAT3912	ST5	5	Stainless steel	187 x 335 x 215	150 x 150 x 300	No	£245.00
BAT3904	ST12	12	Stainless steel	360 x 332 x 215	150 x 300 x 325	Yes	£435.00
BAT3906	ST18	18	Stainless steel	340 x 545 x 215	150 x 300 x 505	Yes	£481.00
BAT3908	ST26	26	Stainless steel	340 x 535 x 270	200 x 300 x 505	Yes	£541.00
BAT3910	ST38	38	Stainless steel	338 x 733 x 260	200 x 300 x 690	Yes	£944.00
BAT3840	P5	5	Plastic	220 x 323 x 180	155 x 160 x 240	No	£170.00
BAT3842	P12	12	Plastic	340 x 412 x 180	155 x 160 x 240	No	£176.00
BAT3844	P18	18	Plastic	340 x 589 x 180	155 x 290 x 510	No	£198.00

Optima Lids

Lids help reduce evaporation/heat loss and avoid sample contamination. Recommended when operating at 60°C and above.



Code	Alt Ref	Description	For use with	Price
BAT3850	PL5	Flat stainless steel	P5 bath	£54.00
BAT3852	PL12	Curved plastic	P12 bath	£46.00
BAT3854	PL18	Curved plastic	P18 bath	£53.00
BAT3920	STL5	Flat stainless steel	ST5 bath	£63.00
BAT3914	STL12	Gabled, hinged (removable) stainless steel	ST12 bath	£223.00
BAT3916	STL26	Gabled, hinged (removable) stainless steel	ST18/ST26 bath	£273.00
BAT3918	STL38	Gabled, hinged (removable) stainless steel	ST38 bath	£408.00