



## Heathrow Scientific Pipette Racks



- Angled four shelved compartments keep your 0.1 to 50mL pipettes or cans within easy reach while preventing damage and clutter
- Heavy duty acrylic version provides crystal-clear viewing from any angle
- ABS versions have cut-out viewing window and come with two removable attachments that hold pipette bulbs, fillers, and pumps
- Magnetised ABS version has strong magnets on both sides and mounts securely to most metal surfaces
- Dims, w x d x h: 88 x 291 x 401mm

| Code           | Alt Ref  | Material       | Colour | Price          |
|----------------|----------|----------------|--------|----------------|
| <b>PIP5510</b> | HS20611  | Acrylic        | Clear  | <b>£154.00</b> |
| <b>PIP5512</b> | HS20612  | ABS            | Blue   | <b>£97.00</b>  |
| <b>PIP5514</b> | HS20613  | ABS            | Red    | <b>£97.00</b>  |
| <b>PIP5518</b> | HS20612M | ABS magnetised | Blue   | <b>£97.00</b>  |



## Heathrow Scientific Universal Pipette Carousels



- Holds up to 8 single channel pipettes, 4 multi channel pipettes, or a combination of both
- Inserts on the upper holders can be added or removed to accommodate different pipette widths
- Rotating head and stable weighted base for one handed removal of pipettes
- Ideal for labs with multiple brands of pipettes
- Conveniently rotate stand head to select or store pipettes
- Quick grip handle to move pipette stand (loaded or unloaded)
- Can be cleaned with soap/water or common laboratory disinfectants (not autoclavable)
- Made of ABS
- Dims, base Ø x top Ø x h: 158 x 229 x 408mm

| Code           | Alt Ref  | Colour        | Price          |
|----------------|----------|---------------|----------------|
| <b>PIP3920</b> | HS120480 | Blue/green    | <b>£122.00</b> |
| <b>PIP3922</b> | HS120481 | Purple/clear  | <b>£127.00</b> |
| <b>PIP3924</b> | HS120499 | Clear/rainbow | <b>£127.00</b> |



## Heathrow Scientific Pipette Tip Bins



- Functional acrylic storage containers for tips, tubes and syringes
- 3mm thick clear acrylic
- Flip-top covers and generous openings allow for easy filling and dispensing
- Non-skid rubber feet

| Code           | Alt Ref | Description   | Dims, l x w x d, mm | Price          |
|----------------|---------|---------------|---------------------|----------------|
| <b>ST08620</b> | HS23401 | 1 compartment | 184 x 194 x 232     | <b>£82.00</b>  |
| <b>PIP9400</b> | HS23402 | 2 compartment | 150 x 177 x 368     | <b>£158.00</b> |
| <b>PIP9402</b> | HS23403 | 3 compartment | 217 x 170 x 217     | <b>£180.00</b> |
| <b>PIP9404</b> | HS23404 | 4 compartment | 287 x 170 x 217     | <b>£145.00</b> |