

Heathrow Scientific Pipette Stands



- Designed to hold most major pipette brands
- Upper slots are 25mm wide and the lower slots are 20mm wide
- Acrylic stands are clear; ABS stands are solid opaque colours

Code	Alt Ref	Capacity	Colour	Material	Price
PIP9452	HS206203	3 pipettes	Clear	Acrylic	£55.00
PIP9454	HS206204	4 pipettes	Clear	Acrylic	£58.00
PIP9456	HS20620	6 pipettes	Clear	Acrylic	£74.00
PIP9460	HS20613C	6 pipettes	Blue	ABS	£84.00



Heathrow Scientific Universal Pipette Wall Mounts



- Can securely hold virtually every known brand of pipette
- Inserts on the upper holders accommodate different pipette widths
- Different coloured upper and lower segments enables colour coding
- Made of ABS
- Easy one handed removal or replacement with the segments unique design that includes traction grips that help centre and hold the pipettes
- Pipettes held in a vertical position to help prevent contamination
- Comes with 2 x wall mounting brackets that can hold 2 pipettes each, and upper and lower segments that easily slide into place



NEW

Code	Alt Ref	Capacity	Colour	Price
PIP0164	HS120575	4 pipettes	Lower holders: 2 x blue, 2 x green, 2 x clear; upper holders: 2 x blue, 2 x green, 1 x purple, 1 x yellow	£37.50

Kartell Automatic Pipette/Burette Rinsing Set



- High speed system taking pipettes up to 600mm long
- Recommended for use where radioactive, caustic and toxic materials are in use
- Constructed of polyethylene which is mounted to protect tips of pipettes from breakage



Code	Alt Ref	Description	Ø, mm	Height, mm	Pack	Price
WAS3100	215	Burette rinser	170	980	Each	£186.00
WAS3102	217	Pipette rinser	170	734	Each	£178.00
WAS3104	218	Pipette jar	162	650	Each	£111.00
WAS3106	219	Basket	145	280	Each	£50.00
WAS3108	1219	Basket (WAS3106) extension handle	—	—	2	£26.00