

Gilson PIPETMAN® Classic Pipettes



The PIPETMAN® Classic line, the world standard in pipetting, offers a wide range of continuously adjustable, air-displacement pipettes with legendary accuracy and precision. Made of highly resistant materials like PVDF and stainless steel, the PIPETMAN® Classic offers uncontested durability and reliability. Built-to-last, it is the economical solution for years of pipetting. Routine cleaning and maintenance are all that is required to keep the PIPETMAN® Classic in top condition for years. The PIPETMAN® Classic design allows for easy access to user-serviceable parts, along with a fully autoclavable tip holder and ejector.

- Gilson's legendary performance
- Robust and built to last
- Easiest pipette to maintain
- Upgraded tip ejector for convenient access to user serviceable parts

PIPETMAN® Classic starter kit includes:

- 3 x PIPETMAN® Classic pipettes (P20, P200, P1000)
- 2 x racks DIAMOND tips (D200, D1000)
- 3 x SINGLE™ pipette hangers
- 'Gilson Guide to Pipetting' handbook

PIPETMAN® Classic microvolume kit includes:


- 3 x PIPETMAN® Classic pipettes (P2, P10, P100)
- 2 x racks DIAMOND tips (DL10, D200)
- 3 x SINGLE™ pipette hangers
- 'Gilson Guide to Pipetting' handbook

PIPETMAN® Classic large volume kit includes:

- 2 x PIPETMAN® Classic pipettes (P5000, P10mL)
- 2 x bags of 10 filters
- 2 x safety bags

PIPETMAN® Classic four pipette kit includes:

- 4 x PIPETMAN® Classic pipettes (P2, P20, P200, P1000)
- 2 x racks DIAMOND tips (D10L, D200, D1000)
- 4 x SINGLE™ pipette holders
- 4 x PIPETMAN® comfort handles
- 'Gilson Guide to Pipetting' handbook



Code	Description	Volume Range, µL	Price
F144801	P2	0.2 - 2	£226.00
F144802	P10	1 - 10	£226.00
F123600	P20	2 - 20	£226.00
F123615	P100	20 - 100	£226.00
F123601	P200	5 - 200	£226.00
F123602	P1000	200 - 1000	£226.00
F123603	P5000	1000 - 5000	£226.00
F161201	P10mL	1000 - 10000	£226.00
F167300	Starter kit	—	£613.00
F167400	Large volume kit	—	£409.00
F167500	Microvolume kit	—	£613.00
F167380	Four pipette kit	—	£794.00