



Gilson PIPETMAN® L Pipettes



A lightweight pipette with a handle that fits naturally in the hand. The PIPETMAN® L features a patented volume locking system for safe, secure pipetting. These features, in combination with Gilson's latest lubricated piston assembly mean that the PIPETMAN® L offers very light pipetting forces.

- Adjustable tip ejector button to adapt to both right- and left-handed users
- Clip mechanism allows easy and safe removal of tip ejector
- GLP features include 2-D barcode for traceability
- Four multichannel models are equipped with V-rings at the bottom of the tip holders for a proper fit with most tip brands

PIPETMAN® L starter kit includes:

- 3 x PIPETMAN® L pipettes (P20L, P200L, P1000L)
- 2 x racks DIAMOND tips (D200, D1000)
- 3 x SINGLE™ pipette holders
- 3 x plastic tip ejectors
- 'Gilson Guide to Pipetting' handbook
- 'Two-minute Inspection' poster

PIPETMAN® L four pipette kit includes:

- 4 x PIPETMAN® L pipettes (P2L, P20L, P200L, P1000L)
- 3 x racks DIAMOND tips (D10L, D200, D1000)
- 4 x SINGLE™ pipette holders
- 'Gilson Guide to Pipetting' handbook

PIPETMAN® L microvolume kit includes:

- 3 x PIPETMAN® L pipettes (P2L, P10L, P100L)
- 2 x racks DIAMOND tips (DL10, D200)
- 3 x SINGLE™ pipette holders
- 3 x plastic tip ejectors
- 'Gilson Guide to Pipetting' handbook
- 'Two-minute Inspection' poster

| Code | Description | Channel | Volume Range, µL | Ejector | Price |
|----------|---------------------------|---------|------------------|---------|---------|
| FA10001M | P2L | Single | 0.2 - 2 | Metal | £249.00 |
| FA10002M | P10L | Single | 1 - 10 | Metal | £249.00 |
| FA10003M | P20L | Single | 2 - 20 | Metal | £249.00 |
| FA10004M | P100L | Single | 10 - 100 | Metal | £249.00 |
| FA10005M | P200L | Single | 20 - 200 | Metal | £249.00 |
| FA10006M | P1000L | Single | 100 - 1000 | Metal | £249.00 |
| FA10001P | P2L | Single | 0.2 - 2 | Plastic | £249.00 |
| FA10002P | P10L | Single | 1 - 10 | Plastic | £249.00 |
| FA10003P | P20L | Single | 2 - 20 | Plastic | £249.00 |
| FA10004P | P100L | Single | 10 - 100 | Plastic | £249.00 |
| FA10005P | P200L | Single | 20 - 200 | Plastic | £249.00 |
| FA10006P | P1000L | Single | 100 - 1000 | Plastic | £249.00 |
| FA10007 | P5000L | Single | 500 - 5000 | — | £249.00 |
| FA10008 | P10mL | Single | 1000 - 10000 | — | £249.00 |
| FA10013 | L 8 channel | 8 | 0.5 - 10 | — | £713.00 |
| FA10009 | L 8 channel | 8 | 2 - 20 | — | £713.00 |
| FA10011 | L 8 channel | 8 | 20 - 200 | — | £713.00 |
| FA10035 | L 8 channel with V-rings | 8 | 20 - 200 | — | £748.00 |
| FA10015 | L 8 channel | 8 | 20 - 300 | — | £713.00 |
| FA10037 | L 8 channel with V-rings | 8 | 20 - 300 | — | £748.00 |
| FA10014 | L 12 channel | 12 | 0.5 - 10 | — | £831.00 |
| FA10010 | L 12 channel | 12 | 2 - 20 | — | £831.00 |
| FA10012 | L 12 channel | 12 | 20 - 200 | — | £831.00 |
| FA10036 | L 12 channel with V-rings | 12 | 20 - 200 | — | £873.00 |
| FA10016 | L 12 channel | 12 | 20 - 300 | — | £831.00 |
| FA10038 | L 12 channel with V-rings | 12 | 20 - 300 | — | £873.00 |
| F167550 | Microvolume kit | Single | — | — | £690.00 |
| F167350 | Starter kit | Single | — | — | £690.00 |
| F167370 | Four pipette kit | Single | — | — | £897.00 |