

Sartorius Proline® Pipettes

SARTORIUS

Manufactured from the latest materials ensuring a light weight and durable pipette and covers the volume range 0.1 to 5000µL.

- Tip ejector that allows convenient one-handed operation
- Finger support keeps the instrument in place with minimum user effort
- Digital display clearly reads volume setting
- Ejector collar and tip cone can be removed for easy cleaning and maintenance
- Durable tip cone materials provide excellent chemical resistance
- Unique safe-cone filter prevents contamination
- Lightweight ergonomic design
- Easy calibration and maintenance
- Colour coded ID stickers included with each pipette



Code	Alt Ref	Channel	Volume Range, µL	Price
PIP7600	720015	Single	0.5 - 10	£197.00
PIP7614	720080	Single	2 - 20	£197.00
PIP7602	720025	Single	5 - 50	£197.00
PIP7608	720050	Single	10 - 100	£197.00
PIP7612	720070	Single	20 - 200	£182.00
PIP7610	720060	Single	100 - 1000	£182.00
PIP7616	720110	Single	1000 - 5000	£182.00
PIP7624	720210	8	0.5 - 10	£574.00
PIP7626	720220	8	5 - 50	£574.00
PIP7630	720240	8	50 - 300	£574.00
PIP7636	720310	12	0.5 - 10	£628.00
PIP7638	720320	12	5 - 50	£628.00
PIP7642	720340	12	50 - 300	£628.00

Accessories

Code	Alt Ref	Description	Price
PIP7675	LH-725630	Carousel stand for 6 pipettes	£101.00
PIP7678	725620	Linear stand for up to 12 pipettes	£69.00

Sustainable laboratory solutions!

We've planted 2500 trees to offset the carbon footprint of this catalogue by double the amount of CO₂, making it climate positive.

Look out for the following logos on our literature:

