

Nickel Electro Clifton® Lids for Stirred Water Baths



Lids help to reduce water evaporation, assists in keeping samples and water free from contamination and saves energy. Stainless steel gabled lids have a large insulated lifting handle.

Code	Alt Ref	For use with	Price
BAT7095	SL4-8	NE4-8P, NE4-8D, NE4-8T, NE4-8HT	£140.00
BAT7098	SL4-14	NE4-14P, NE4-14HT, NE4-14D, NE4-14T	£144.00
BAT7100	SL4-22	NE4-28D, P04, NE4-28P, NE4-22HT, NE4-28T, NE4-22P, NE4-22T, NE4-22D	£162.00
BAT7162	SL1-22H	NE5-28D, NE1D-22, NE2-28D, NE1-28, NE3-28DT, NE5-28, NE2-22D, NE1-22	£273.00



Grant Aqua Pro Shaking Water Baths



Offers exceptional temperature performance and usability. With three bath models to choose from, there is a shaking water bath to meet basic and demanding applications in research, industry and life sciences. Includes UK and EU cord sets and clear polycarbonate lid. Universal tray included as standard, other shaking trays sold separately.

- Temperature range: ambient +5 to 99°C
- Stability: $\pm 0.1^{\circ}\text{C}$
- Digital control/LED display
- Linear shaking speed 20 – 200RPM
- Adjustable shaking speed and stroke length (BAT1056 model OLS26 only)
- Countdown timer (1 – 999 mins)
- Soft start and extremely quiet during operation
- Front located drain tap
- User defined high temperature alarm and cut off



Code	Model	Shaking Mode	Capacity, L	Ext Dims, w x d x h, mm	Tray Area, w x d, mm	Price
BAT1056	OLS26	Orbital & linear	26	335 x 590 x 475	380 x 235	£2650.00
BAT1058	LSB12	Linear	12	360 x 385 x 425	240 x 235	£1890.00
BAT1060	LSB18	Linear	18	335 x 565 x 425	420 x 235	£2150.00

Accessories

Code	Alt Ref	Description	For Model	Price
BAT1090	CC26	Coiling coil for bath operation operated at or below ambient, down to 0°C	OLS26	£1850.00
BAT1092	CW26	Heat exchange coil attach to water supply or refrigerated circulator, can be used down to 2°C above the temp of the coolant	OLS26	£220.00

View our chemicals catalogue

Visit the SLS website to read online or request your copy:

scientificlabs.co.uk/literature (UK & Northern Ireland)
scientificlabs.ie/literature (Republic of Ireland)
scientificlabs.com/literature (International)

