

Sartorius Proline® Plus Pipettes

SARTORIUS

Offers comfort and quality for your everyday manual pipetting. It combines durable construction with ease and lightness of use. Easy volume adjustment with click stop system and big, forward facing, volume display help the user always to choose correct volume.

- Low pipetting forces, a comfortable handle and an ergonomic finger support for effortless pipetting
- Contamination free pipetting with replaceable Safe-Cone Filters
- Spring loaded tip cones with Optiload® mechanism which allows tips to be loaded and ejected with minimum force
- No opening tools are needed for cleaning and maintenance
- Robust design and easy to maintain
- Fully autoclavable



Code	Alt Ref	Channel	Volume Range, µL	Price
PIP3500	728010	Single	0.1 - 3	£206.00
PIP3502	728020	Single	0.5 - 10	£206.00
PIP3504	728030	Single	2 - 20	£206.00
PIP3506	728050	Single	10 - 100	£206.00
PIP3508	728060	Single	20 - 200	£206.00
PIP3510	728070	Single	100 - 1000	£206.00
PIP3512	728080	Single	500 - 5000	£206.00
PIP3514	728090	Single	1000 - 10000	£206.00
PIP3516	728120	8	0.5 - 10	£601.00
PIP3518	728130	8	10 - 100	£601.00
PIP3520	728140	8	30 - 300	£601.00
PIP3522	728220	12	0.5 - 10	£655.00
PIP3524	728230	12	10 - 100	£655.00
PIP3526	728240	12	30 - 300	£655.00

Multipacks

Code	Alt Ref	Description	Contains	Price
PIP1548	LH-728670	Proline Plus 3 -pack 10	3 x single channel pipettes [0.5-10/10-100/100-1000µL]; 3 x tip racks - 96 tips [0.1-10/0.5-200/10-1000µL]; linear stand; 3 x pipette holders	£493.00
PIP1550	LH-728671	Proline Plus 3 -pack 20	3 x single channel pipettes [2-20/20-200/100-1000µL]; 2 x tip racks - 96 tips [0.5-200/10-1000µL]; linear stand; 3 x pipette holders	£493.00
PIP1552	LH-728672	Proline Plus 4 -pack 20	4 x single channel pipettes [0.5-10/2-20/20-200/100-1000µL]; 3 x tip racks - 96 tips [0.1-10/0.5-200/10-1000µL]; linear stand; 4 x pipette holders	£643.00
PIP1554	LH-728673	Proline Plus 4-pack 100	4 x single channel pipettes [0.5-10/10-100/20-200/100-1000µL]; 3 x tip racks - 96 tips [0.1-10/0.5-200/10-1000µL]; linear stand; 4 x pipette holders	£643.00
PIP1556	LH-728674	Proline Plus 5-pack 10	5 x single channel pipettes [0.5-10/10-100/20-200/100-1000/500-5000µL]; 2 x tip racks - 96 tips [0.5-200/10-1000µL]; 1 x tip rack - 100 tips [5000µL]; linear stand; 5 x pipette holders	£797.00
PIP1558	LH-728675	Proline Plus 5 -pack 20	5 x single channel pipettes [2-20/10-100/20-200/100-1000/500-5000µL]; 2 x tip racks - 96 tips [0.5-200/10-1000µL]; 1 x tip rack - 100 tips [5000µL]; linear stand; 5 x pipette holders	£797.00
PIP1560	LH-728676	Proline Plus 3+1 -pack	3 x single channel pipettes [0.5-10/20-200/100-1000µL]; 1 x 8 channel pipette [30-300µL]; 4 x tip racks - 96 tips [0.1-10/0.5-200/5-350/10-1000µL]; linear stand; 4 x pipette holders	£892.00
PIP1562	LH-728677	Proline Plus 3 -pack Max	3 x single channel pipettes [100-1000/500-5000/1000- 10000µL]; 1 x tip racks - 96 tips [10-1000µL]; 2 x Bulk in a Bag tips [100 x 0.1-5mL/250 x 1-10mL]; linear stand; 3 x pipette holders	£493.00