

## Eppendorf Reference® 2 Pipettes

eppendorf

Reference® 2 is the successor to the Reference® series of pipettes incorporating a new design, reduced weight and operating forces and now including multi-channel versions. Reference® 2 provides the highest levels of accuracy with reliable handling and user safety making Reference® 2 ideal for use with precious liquids, or in any application requiring high levels of accuracy.


- One button operation ensure fast and ergonomic handling, reduced operating effort and active aerosol reduction
- Spring loaded tip cone improves ergonomics by lowering tip attachment force
- The Reference® 2 multi channel includes a spring-loaded tip cone which can be switched on/off
- Quick lock enables easy cleaning and maintenance of your pipette
- 4 digit volume display with magnifying window for easier volume identification
- Robust with external edges made from stainless steel
- Embedded RFID chip contains all relevant data regarding the pipette

Code	Channel	Volume, µL	Price
E4924000010	Single	0.1 - 2.5	£313.00
E4924000029	Single	0.5 - 10	£313.00
E4924000037	Single	2 - 20	£313.00
E4924000053	Single	10 - 100	£313.00
E4924000061	Single	20 - 200	£313.00
E4924000070	Single	30 - 300	£313.00
E4924000088	Single	100 - 1000	£313.00
E4924000096	Single	250 - 2500	£313.00
E4924000100	Single	500 - 5000	£314.00
E4924000118	Single	1000 - 10000	£314.00
E4926000018	8	0.5 - 10	£862.00
E4926000026	12	0.5 - 10	£944.00
E4926000034	8	10 - 100	£862.00
E4926000042	12	10 - 100	£944.00
E4926000050	8	30 - 300	£862.00
E4926000069	12	30 - 300	£944.00



### Starter Packs

Code	Description	Supplied with	Price
E4924000908	Starter pack option 1	epT.I.P.S. and 0.5-10µL, 10-100µL, 100-1000µL pipettes	£733.00
E4924000916	Starter pack option 2	epT.I.P.S. and 2-20µL, 20-200µL, 100-1000µL pipettes	£733.00
E4924000924	Starter pack option 3	epT.I.P.S. and 100-1000µL, 0.5-5mL, 1-10mL pipettes	£733.00



## Eppendorf Research® plus Pipettors

eppendorf

Eppendorf sets a new standard in pipetting. The ultra-light Eppendorf Research® plus pipette meets the highest needs in precision and accuracy - combined with ultimate ergonomics and increased flexibility.

- Light weight for ease of use and comfort
- Low tip-ejection forces for stress-free workflow
- Fully autoclavable to save time