



## Corning® Lambda™ EliteTouch™ Pipette Starter Kit

CORNING

Engineered to provide the highest level of ergonomic, accuracy, and precision. Lightweight construction, contoured shape of the handle, and 4-digit counter were designed to guarantee comfortable pipetting. Smooth plunger movement and extremely low pipetting forces reduce wrist strain and fatigue (RSI). Coloured pushbuttons help define user or application for sample safety and lower risk of cross contamination.

- Ergonomic handle design for perfect fit for both right and left hand
- Convenient one-hand volume setting with auto-lock system
- 4 digit counter for greater precision, perfectly positioned for visibility at all times
- Fully autoclavable and UV resistant
- Easy in-lab calibration helps you quickly adjust your pipettor to various liquid types
- Individually tested and supplied with a certificate of quality

### Kit includes:

- 4 x Lambda™ EliteTouch™ pipettes (0.5 - 10µL, 2 - 20µL, 20 - 200µL, 100 - 1000µL)
- 1 x universal linear stand for 4 single channel pipettes
- 3 x DeckWorks™ pipette tips (10µL, 200µL, 1000µL)
- 4 x 3 coloured pushbuttons

Code	Price
6065	£928.00



## Sartorius Picus® Electronic Pipettes

SARTORIUS

Picus® offers all the advantages you would expect from a Sartorius pipette whilst the NxT version is targeted at customers working in a regulated environment where full compliance is a necessity.

### Picus® Electronic Pipettes

Incorporating an electronic brake and piston control system help ensure accurate and precise pipetting results whilst an intuitive user interface with an adjustment wheel ensure fast volume setting and easy menu navigation.

- Electronic tip ejection
- Built in plate tracker functionality guides pipetting steps for improved microplate work
- Light at 100g
- Optiload ensures tip loading with a constant and equal force for perfect sealing
- Hot key access for up to 10 saved programmes

Code	Alt Ref	Channel	Volume Range, µL	Price
PIP3800	735021	Single	0.2 - 10	£624.00
PIP3802	735041	Single	5 - 120	£624.00
PIP3804	735061	Single	10 - 300	£624.00
PIP3806	735081	Single	50 - 1000	£624.00
PIP3808	735101	Single	100 - 5000	£624.00
PIP3810	735111	Single	500 - 10000	£624.00
PIP3812	735321	8	0.2 - 10	£965.00
PIP3814	735341	8	5 - 120	£965.00
PIP3816	735361	8	10 - 300	£965.00
PIP3818	735391	8	50 - 1200	£965.00
PIP3820	735421	12	0.2 - 10	£1151.00
PIP3822	735441	12	5 - 120	£1151.00
PIP3824	735461	12	10 - 300	£1151.00
PIP3826	735491	12	50 - 1200	£1151.00