

Eppendorf Xplorer® and Xplorer® plus
Pipettes Continued

eppendorf

Accessories

Code	Description	For use with	Price
E3116000023	Charger Carousel 2	6 x Xplorer® or Xplorer® plus	£410.00
E3116000031	Charger Stand 2	1 x Xplorer® or Xplorer® plus	£217.00

Eppendorf Xplorer® plus Move It® Pipettors
with Adjustable Tip Spacing

eppendorf

Do you frequently transfer samples between different vessel formats, such as between tubes and plates, and strive for more efficient and safe solutions? The new adjustable tip spacing pipette Move It® accelerates and simplifies synchronous pipetting of a series of samples. A unique design without tubing connections results in increased precision and durability as well as user and sample safety.

- Efficient sample transfers from one vessel type to another.
- Pre-defined settings allow for quick switches back and forth between different formats
- Tip spacing freely selectable between 4.5 and 33mm
- Time-saving thanks to rotatable pipette head – 360°
- Tubeless system for reliability and reproducible results- no cable tangle and tubing-porosity, for consistently reliable robustness and precision
- Charger included



Code	Channel	Volume Range, µL	Tip Spacing Range, mm	Price
E4861000816	4	15 - 300	9 - 33	£1596.00
E4861000833	4	50 - 1200	9 - 33	£1596.00
E4861000817	6	15 - 300	9 - 20	£1646.00
E4861000834	6	50 - 1200	9 - 20	£1646.00
E4861000781	8	1 - 20	4.5 - 14	£1696.00
E4861000794	8	5 - 100	4.5 - 14	£1696.00
E4861000818	8	15 - 300	9 - 14	£1696.00
E4861000835	8	50 - 1200	9 - 14	£1696.00
E4861000782	12	1 - 20	4.5 - 9	£1796.00
E4861000795	12	5 - 100	4.5 - 9	£1796.00

Eppendorf VisioNize® Pipette Manager

eppendorf

- Manage your pipette fleet more effectively and enhance collaboration in teams. Integrated features like Pipette Set Creation or Liquid Type Sharing support smooth and efficient operations
- Easily convert your conventional Eppendorf Xplorer® and Eppendorf Xplorer® plus electronic pipettes into a connected electronic pipette with the WiFi module
- Quicker input of volumes via the VisioNize® pipette manager, all settings will be transferred to connected electronic pipettes for an immediate start of pipetting
- Get guidance by the software-embedded Liquid Manager when working with liquids of different viscosity, volatility, or surface tension, or foam-building liquids



Code	Description	Price
E1004000001	VisioNize® pipette manager	£839.00
E4861000970	Xplorer Connect® WiFi module inc. battery	£70.00