



## Gilson PIPETMAN® M Connected Pipettes



The PIPETMAN® M Connected Bluetooth®-enabled electronic pipette offers accuracy and precision in both standard and repetitive pipetting modes. The motorised piston allows for high pipetting reproducibility and accuracy in all modes. Requiring virtually no pipetting forces to aspirate and dispense samples, PIPETMAN® M Connected helps reduce pipetting fatigue and increase pipetting efficiency.

Connect the pipette to the Gilson PipetteScope™ app and know exactly when the pipette needs to be calibrated, access service history reports and gain visibility of the pipette usage analytics. With the PIPETMAN® M Connected software, create custom protocol from 11 different pipetting task on a PC or Android and wirelessly sync up to ten custom protocols to your pipette. Used together with TRACKMAN® Connected, the pipette interacts in real time with device to plan, execute and track pipetting protocols faster, easier and more reliable minimising human error.

- Intuitive motorised pipette with Bluetooth® connectivity
- Accurate and precise pipetting
- On-target performance
- Access pipettes information via Gilson PipetteScope™ app

### PIPETMAN® M single channel starter kit includes:

- TRACKMAN® Connected tablet and accessories
- 3 x PIPETMAN® M Connected pipettes (P200M, P1200M, P5000M)
- 1 x rack of tips (5mL)
- 10 x racks of filter tips (5 x 500µL, 5 x 1200µL)

### PIPETMAN® M starter kit pro includes:

- TRACKMAN® Connected tablet and accessories
- 3 x PIPETMAN® M Connected pipettes (P200M, P1200M, P5000M, P8x1200M)
- 1 x rack of tips (5mL)
- 15 x racks of filter tips (5 x 200µL, 10 x 1200µL)

### PIPETMAN® M multi channel starter kit includes:

- TRACKMAN® Connected tablet and accessories
- 3 x PIPETMAN® M Connected pipettes (P1200M, P8x200M, P8x1200M)
- 20 x racks of filter tips (10 x 200µL, 10 x 1200µL)

Code	Description	Channel	Volume Range, µL	Price
F81040	P10M	Single	0.5 - 10	£498.00
F81041	P20M	Single	2 - 20	£498.00
F81042	P100M	Single	5 - 100	£498.00
F81043	P200M	Single	20 - 200	£498.00
F81044	P300M	Single	20 - 300	£498.00
F81045	P1200M	Single	100 - 1200	£498.00
F81046	P5000M	Single	500 - 5000	£498.00
F81047	P10mLM	Single	1 - 10000	£498.00
F81048	P8x10M	8	0.5 - 10	£918.00
F81050	P8x20M	8	1 - 20	£918.00
F81052	P8x100M	8	10 - 100	£918.00
F81054	P8x200M	8	20 - 200	£918.00
F81056	P8x300M	8	10 - 300	£918.00
F81058	P8x1200M	8	50 - 1200	£918.00
F81049	P12x10M	12	0.5 - 10	£1031.00
F81051	P12x20M	12	1 - 20	£1031.00
F81053	P12x100M	12	10 - 100	£1031.00
F81055	P12x200M	12	20 - 200	£1031.00
F81057	P12x300M	12	10 - 300	£1031.00
F81059	P12x1200M	12	50 - 1200	£1031.00
FB10201	Single starter kit	—	—	£2582.00
FB10211	Multi starter kit	—	—	£3362.00
FB10221	Starter kit pro	—	—	£3401.00