

Gilson TRACKMAN® Connected



A tablet with exclusive accessories and applications that make pipetting on 96 and 384 well plates faster and more reliable, improving efficiency at the bench by tracking pipetting tasks. Designed to communicate with Gilson's PIPETMAN® M Connected via Bluetooth®, the tablet interacts in real time with the pipette and guides the user through their protocol with the PipettePilot™ app. The app displays where, when, and how much to pipette and even captures notes, pictures and monitors the protocol. It can even adjust the pipette volume automatically according to the method.


- Improves traceability and reproducibility
- Accelerates collaborative reporting
- Reports on your pipetting performance
- Guides you through your pipetting protocol
- Increases your efficiency and reduces costly errors



Code	Price
FB1021	£1586.00

Accessories

Code	Description	Price
FB02201	Tablet holder	£102.00
FB02202	Plate holder	£26.00



Gilson MICROMAN® E Pipettes



The ergonomically designed MICROMAN® E positive displacement pipettes, along with disposable capillary piston (CP) tips, provide the highest precision when pipetting viscous, non-aqueous liquids such as oils, blood, glycerol, and more, while protecting the pipette from contamination.

- Relies on positive-displacement technology to improve the precision of transferring viscous and volatile liquid
- Improved handle shape that rests comfortably in your hand, with a large push button for easy aspiration and dispensing
- Volume-locking switch allows for easy volume adjustment and minimises the risk of accidentally changing the volume
- Ideal when dealing with infectious solutions and radiolabeled components



Code	Description	Volume Range, µL	Price
FD10001	M10E	1 - 10	£323.00
FD10002	M25E	3 - 25	£323.00
FD10003	M50E	20 - 50	£323.00
FD10004	M100E	10 - 100	£323.00
FD10005	M250E	50 - 250	£323.00
FD10006	M1000E	100 - 1000	£323.00