

## Nickel Electro Clifton® NE5 Shaking Water Baths



Available with a broad range of accessory trays, flask clips and racks to provide flexibility for routine and specialised work. 10L bath includes perforated accessory tray. Accessory trays for 28L baths are sold separately.


- Choice of analogue or PID digital control
- Temperature range ambient +5 to 99°C
- Easy to use splash proof controls
- Linear shaker
- Features auto setting over and under temperature alarms
- Accessory clips and racks can be mounted onto perforated shelf to accommodate flasks, micro titre plates or test tube racks
- High quality corrosion resistant stainless steel tank housed in a durable outer case
- Painted surfaces feature antibacterial finish
- Heaters are concealed under the tank providing a clutter free working area easy to wipe clean
- Every water bath undergoes a factory calibration at a reference temperature validating each bath for total confidence in use



| Code           | Alt Ref | Capacity, L | Control  | Ext Dims, h x d x w, mm | Platform Dims, mm | Shaking Speed, RPM | Price           |
|----------------|---------|-------------|----------|-------------------------|-------------------|--------------------|-----------------|
| <b>BAT7148</b> | NE5-10D | 10          | Digital  | 290 x 449 x 332         | 175 x 270         | 0 - 400            | <b>£1988.00</b> |
| <b>BAT7149</b> | NE5-28  | 28          | Analogue | 290 x 715 x 332         | 376 x 276         | 0 - 380            | <b>£2061.00</b> |
| <b>BAT7150</b> | NE5-28D | 28          | Digital  | 290 x 715 x 332         | 376 x 276         | 0 - 400            | <b>£2181.00</b> |

### Lids

Lids help to reduce water evaporation, assists in keeping samples and water free from contamination and saves energy. Clear polycarbonate gabled lids can be hinged back out of the way to access the water bath. Stainless steel gabled lids have a large insulated lifting handle.



| Code           | Alt Ref | Description                     | For use with   | Price          |
|----------------|---------|---------------------------------|--|----------------|
| <b>BAT7071</b> | GL1-8   | Hinged polycarbonate gabled lid | NE2-10D, NE2-8D, NE1-8, NE5-10D  | <b>£73.00</b>  |
| <b>BAT7094</b> | SL1-22  | Stainless steel gabled lid      | NE5-28D, NE2-28D, NE1-9, NE1D-22, NE1-28, NE2D-8/22, NE2-9D, NE3-28DT, NE5-28, NE2-22D, NE1-22 | <b>£159.00</b> |
| <b>BAT7162</b> | SL1-22H | Stainless steel gabled lid      | NE5-28D, NE1D-22, NE2-28D, NE1-28, NE3-28DT, NE5-28, NE2-22D, NE1-22                           | <b>£273.00</b> |
| <b>BAT7160</b> | GL1-22  | Hinged polycarbonate gabled lid | NE2-28D, NE5-28, NE1-22, NE2-22D, NE1-28, NE1B-28, NE5-28D, NE1D-22                            | <b>£73.00</b>  |