

Gilson PIPETMAN® DIAMOND Pipette Tips Continued



TOWERPACK™ Tips

- A quick, easy and affordable rack refill system
- Uses 60% less space on your bench and 70% less plastic than regular rack box
- Whole range of TOWERPACK™ fits the amazingly easy to use universal reload box
- One hand system operation
- Fully autoclavable

Code	Description	Volume Range, µL	Sterility	Pack	Price
F167101	D10	0.1 - 10	Non sterile	10 x 96	£39.00
F167102	DL10	0.1 - 20	Non sterile	10 x 96	£39.00
F167103	D200	2 - 200	Non sterile	10 x 96	£34.00
F167104	D1000	100 - 1000	Non sterile	7 x 96	£27.00
F167105	D300	20 - 300	Non sterile	10 x 96	£34.00
F167201	D10	0.1 - 10	Sterile	10 x 96	£49.00
F167202	DL10	0.1 - 20	Sterile	10 x 96	£49.00
F167203	D200	2 - 200	Sterile	10 x 96	£44.00
F167204	D1000	100 - 1000	Sterile	7 x 96	£34.00
F167205	D300	20 - 300	Sterile	10 x 96	£44.00



A

Reload Pack Tips

- Quickest way to reload empty TIPACK™ box
- Uses 70% less plastic, recyclable packaging is eco-friendly
- Three times smaller than a regular tips box, it is easy to store on your bench
- Affordable and convenient
- Racks are fully autoclavable

Code	Description	Volume Range, µL	Sterility	Pack	Price
F167012	DL10	0.1 - 20	Non sterile	10 x 96	£35.00
F167013	D200	2 - 200	Non sterile	10 x 96	£31.00
F167014	D1000	10 - 1000	Non sterile	6 x 96	£21.00



A

Blister Refill Pack Tips

- Individually wrapped racks of tips that are convenient, clean and eco-friendly
- Ready to use, just peel the lid and insert into the empty box
- Uses 45% less plastic than regular rack box
- Racks are fully autoclavable
- Colour coded for easy identification and compliance with Good Laboratory Practice (GLP)
- Available as autoclavable, sterilised, or sterilised filter tips

Code	Description	Volume Range, µL	Sterility	Pack	Price
F172200	DL10	0.1 - 20	Non sterile	10 x 96	£43.00
F172201	DL10	0.1 - 20	Sterile	10 x 96	£51.00
F172300	DL200	2 - 200	Non sterile	10 x 96	£38.00
F172301	DL200	2 - 200	Sterile	10 x 96	£44.00
F172400	D300	20 - 30	Non sterile	10 x 96	£38.00
F172401	D300	20 - 30	Sterile	10 x 96	£44.00
F172500	D1000	100 - 1000	Non sterile	10 x 96	£40.00
F172501	D1000	100 - 1000	Sterile	10 x 96	£49.00
F172800	D1200	50 - 1200	Non sterile	10 x 96	£50.00
F172801	D1200	50 - 1200	Sterile	10 x 96	£56.00



A