

Gilson PIPETMAN® DIAMOND Pipette Tips Continued



ECOPACK™ Tips

- Purchase in bulk and save money by simply reloading these tips into empty racks
- Bags are sealed with rapid-sealing strip for easy access and to prevent contamination
- Minimal waste



A

Code	Description	Volume Range, µL	Sterility	Pack	Price
F161630	D10	0.1 - 10	Non sterile	10 x 1000	£275.00
F161450	DL10	0.1 - 20	Non sterile	10 x 1000	£275.00
F161930	D200	2 - 200	Non sterile	10 x 1000	£273.00
F161730	D300	20 - 300	Non sterile	10 x 1000	£273.00
F161670	D1000	100 - 1000	Non sterile	10 x 1000	£273.00
F161110	D1200	100 - 1200	Non sterile	10 x 1000	£273.00

Empty Racks/Boxes

Code	Description	For use with	Pack	Price
F167100	Universal reload box	TOWERPACK™	Each	£11.00
F171204	Empty box	Reload/blister packs DL10, D200, D300	2	£10.00
F171504	Empty box	Reload/blister packs D1000, D1200	2	£10.00



Gilson AmpliPur® Expert Filter Tips



ISO 8655 compliant and provide a significant value for the money spent, manufactured according to the certified Gilson quality standard. These filter tips protect precious samples and your pipettes from cross-contamination and are designed to fit various pipette brands. AmpliPur® Expert Tips are certified DNase- and RNase-free, making them the preferred tool for molecular biology applications, supporting you in ensuring the limit of detection of your experiment, resulting in more sensitive nucleic acid quantitation.

- Colour-coded for easy identification
- Compliant with Good Laboratory Practice (GLP)
- Designed to fit various pipette brands for maximum convenience on the bench

Supplied in packs of 960 (10 x 96)



Code	Description	Volume Range, µL	Price
F174101	FT10	0.5 - 10	£79.00
F174201	FT20	1 - 20	£73.00
F174301	FT200	10 - 200	£73.00
F174401	FT1000	100 - 1000	£75.00
F174501	FT1200	100 - 1200	£75.00