

JULABO CORIO™ Refrigerated Circulators

Julabo

The CORIO™ CD series of refrigerated circulators are ideal for standard tasks and routine work in laboratories. With a bright display which can be read at a distance and one click control for both circulation and refrigeration, high quality insulation minimises energy loss and a lockable power plug eliminates the risk of accidental disconnection.

- Very quiet
- Easy changeover from internal to external circulation and vice versa
- External pump connections (M16x1)
- USB interface allows remote control
- Cooling coil design yields more usable space in the bath tank
- Bath lid and drain tap included
- Class III (FL) according to DIN 12876-1
- Precise temperature control
- Easy operation
- Temperature stability: $\pm 0.03^{\circ}\text{C}$
- Heating capacity: 2kW
- Pump capacity: 15L/min, pump pressure: 0.35 bar
- 1PLUS warranty offers an extension to the normal warranty of up to 24 months, limited to a maximum of 10000 working hours



Code	Alt Ref	Model	Temp Range, °C	Cooling Capacity (@ 20°C), kW	Dims, w x d x h, mm	Price
BAT8100	9012701	CD-200F	-20 to +150	0.22	230 x 390 x 650	£3779.00
BAT8101	9012702	CD-201F	-20 to +150	0.22	440 x 410 x 440	£3779.00
BAT8102	9012703	CD-300F	-25 to +150	0.31	240 x 420 x 660	£4395.00
BAT8103	9012704	CD-600F	-25 to +150	0.60	330 x 470 x 690	£5525.00
CD-900F	9012706	CD-900F	-40 to +150	0.90	390 x 620 x 750	£5839.00

SLS Lab Pro Analogue Ultrasonic Baths

LAB PRO

With capacities from 1.5L to 25L, the SLS Lab Pro analogue ultrasonic baths offer simple analogue timer and temperature controls, providing entry level ultrasonic baths ideal for cleaning small items for use in a wide range of applications. Supplied with a stainless steel basket and plastic lid. SLS2006 has a non adjustable heater element.

- Ambient to 70°C temperature control on heated models
- 0 – 15 minute timer
- Compact and user friendly



Code	Capacity, L	Heated	Overall Dims, w x d x h, mm	Int Dims, w x d x h, mm	Price
SLS2004	1.5	No	180 x 160 x 190	150 x 140 x 100	£305.00
SLS2006	1.5	Yes	180 x 160 x 190	150 x 140 x 100	£356.00
SLS2016	2.5	No	270 x 170 x 210	240 x 140 x 100	£390.00
SLS2018	2.5	Yes	270 x 170 x 210	240 x 140 x 100	£451.00
SLS2020	4.5	Yes	330 x 180 x 275	300 x 150 x 150	£801.00
SLS2022	13	Yes	375 x 345 x 275	330 x 300 x 150	£1015.00
SLS2026	25	Yes	550 x 345 x 365	500 x 300 x 200	£1367.00