

Corning® Reagent Reservoirs



- Manufactured from modified polystyrene
- White
- Sterile
- Disposable

| Code | Capacity, mL | Inner Pack | Pack | Price |
|------|-------------------------|------------|------|---------|
| 4871 | 50 | 1 | 100 | £279.00 |
| 4870 | 50 | 5 | 200 | £241.00 |
| 4873 | 100 | 1 | 100 | £329.00 |
| 4872 | 100 | 5 | 200 | £288.00 |
| 4877 | 12 channel (approx 5mL) | 1 | 100 | £195.00 |
| 4878 | 8 channel (approx 7mL) | 1 | 100 | £190.00 |



Blue-Ray Biotech EzMate™ Automatic Pipetting Systems



Ideal for routine low volume pipetting. It helps eliminate tedious pipetting and human error, which enables more precise and standardised results to improve productivity. With highly-customised adapters and application-oriented software, EzMate™ can satisfy the various needs of assays, including PCR/qPCR set up, sample normalisation, NGS preparation, and reagent kit manufacturing methods.

- Easy to use system setup software without need of additional technical support
- Built in automatic PCD/qPCR set-up protocols that support USB data transfer
- Single, 4 and 8 channel, 50µL or 200µL automated pipetting modules (APMs)
- 6 interchangeable standard microplate/tip rack adapters and 2 interchangeable reagent vessel adapters
- Active cooling and heating module (ACHM) options for sensitive samples/reagents
- UV lamp and HEPA filter meet the need of high-contamination assays
- Customised adapters support various types of labware
- Excellent results for qPCR standard curve and replicates
- Better precision and pipetting results than manual pipetting
- A wide variety of application-oriented method support, including PCR/qPCR set up, NGS preparation, sample normalisation, reagent kit manufacturing, cell-based assay, etc.
- Supplied with 15.6" notebook computer and EzStarter control software



| Code | Alt Ref | Description | Dims, w x d x h, mm | Price |
|---------|----------------|--|---------------------|-----------|
| PIP1922 | AREZ401-100000 | EzMate 401 4-position | 590 x 450 x 470 | £11242.00 |
| PIP1924 | EZ4SU00-100000 | EzMate 401s 4-position with UV and HEPA | 660 x 450 x 490 | £12795.00 |
| PIP1926 | EZ4SUA0-100000 | EzMate 401s 4-position with UV and HEPA and 2-position active cooling/heating module | 660 x 450 x 490 | £17764.00 |
| PIP1928 | EZ4SUAC-100000 | EzMate 401s 4-position with UV and HEPA and 3-position active cooling/heating module | 660 x 450 x 490 | £19627.00 |
| PIP1930 | AREZ601-100000 | EzMate 601s 6-position | 740 x 450 x 470 | £13416.00 |
| PIP1932 | EZ6SU00-100000 | EzMate 601s 6-position with UV and HEPA | 810 x 450 x 490 | £15280.00 |
| PIP1934 | EZ6SUA0-100000 | EzMate 601s 6-position with UV and HEPA and 2-position active cooling/heating module | 810 x 450 x 490 | £20248.00 |
| PIP1894 | EZ6SUAC-100000 | EzMate 601s 6-position with UV and HEPA and 3-position active cooling/heating module | 810 x 450 x 490 | £22112.00 |

Accessories

Please enquire for full range of accessories.

| Code | Alt Ref | Description | Price |
|---------|----------------|--|----------|
| PIP1896 | AREZ401-100000 | 1 channel 50µL pipetting module (new bracket) | £932.00 |
| PIP1902 | EZ4SU00-100000 | 1 channel 200µL pipetting module (new bracket) | £932.00 |
| PIP1908 | EZ4SUA0-100000 | 4 channel 50µL pipetting module (new bracket) | £931.00 |
| PIP1906 | EZ4SUAC-100000 | 4 channel 200µL pipetting module (new bracket) | £1180.00 |
| PIP1898 | EZ6SUAC-100000 | 8 channel 50µL pipetting module (new bracket) | £1180.00 |
| PIP1900 | AREZ601-100000 | 8 channel 50µL pipetting module for Axygen FX-50-R (new bracket) | £1180.00 |
| PIP1904 | EZ6SU00-100000 | 8 channel 200µL pipetting module (new bracket) | £1180.00 |