



## Joint Clips for Wide Neck Flasks



- For securing flat flange lids to wide neck reaction flasks and culture vessels
- Stainless steel construction

Code	Flange Bore, mm	To Fit Joint Size	Price
<b>QJC75F</b>	75	FG75	£30.70
<b>QJC100F</b>	100	FG100	£30.80



## Joint Clips for Conical Glass Joints

- For retaining glassware with QUICKFIT® conical joints
- Can considerably reduce the need for clamping
- Easy to fit and remove
- Can be used at temperatures up to 80°C
- Clips are colour coded for easy identification of joint sizes

Supplied in packs of 10

Code	Joint Size	Colour Code	Price
<b>QKC100</b>	10/19	Mint	£41.00
<b>QKC120</b>	12/21	Violet	£34.90
<b>QKC140</b>	14/23	Yellow	£41.00
<b>QKC190</b>	19/26	Blue	£39.00
<b>QKC240</b>	24/29	Green	£27.00
<b>QKC290</b>	29/32	Red	£41.00
<b>QKC400</b>	40/38	Ochre	£47.90
<b>QKC450</b>	45/40	Brown	£69.00



## Screwthread Connectors



- Polypropylene screwthread connectors with Viton 'O' rings for use up to 140°C
- For use with glassware having 13mm screwthread inlets and outlets
- Allows 9mm bore rubber hose and flexible tubing to be attached to the glassware easily and safely

Code	Pack	Price
<b>Q451002</b>	10	£60.00