

Distillation Apparatus



- Combines a normally separate Liebig condenser, stillhead and receiver adapter into a one piece assembly
- Angle and position pre-defined
- Simple to set up with fewer clamps required
- Less grease contamination
- No air leakage under vacuum
- No liquid hold-up at the joints


Code	Cone Size	Socket Size	Approx. Effective Condenser Length, mm	Price
QDA122	19/26	14/23	155	£115.00



Kjeldahl Distillation Assembly

- Convenient assembly of QUICKFIT[®] glassware for Kjeldahl digestions and distillations
- Flask can be used alone for digestion stage and then fitted to rest of assembly for the distillation process
- Manufactured from heat resistant and chemically durable borosilicate glass

Code	Description	Socket Size	Cone Size	Pack	Price
Q280MC	Complete assembly	—	—	Each	£421.00
QC112SC	Condenser	19/26	19/26	Each	£38.50
QD312	Dropping funnel, 50mL	—	19/26	5	£427.00
QFK5003L	Kjeldahl flask, 500mL	24/29	—	Each	£28.70
QMA13	Multiple adapter	19/26	24/29	Each	£42.20
QRA122	Delivery adapter	19/26	—	Each	£21.46
QSH712	Splash head	19/26	19/26	Each	£66.00




Still Head Adapters



- Used in distillation assemblies for connecting flasks to condenser
- Upper socket designed to accept thermometer either via cone/screwthread adapter or thermometer pocket

Code	Cone Size to Fit Flask	Cone Size to Fit Condenser	Socket Size	Price
QSH41	14/23	14/23	14/23	£28.30
QSH42	19/26	19/26	19/26	£41.20
QSH433	24/29	24/29	14/23	£35.40



Liquid handling... Down to the last drop!

From pipettes to dispensers, tips to pipetting stations, SLS is your No.1 source for liquid handling.

See pages 745–815



**SCIENTIFIC
LABORATORY
SUPPLIES**

