

## Heathrow Scientific Adapt-a-Rack® Continued



Code	Alt Ref	Colour	Price
RAC0130	HS120184	Blue/green	£16.10
RAC0132	HS120185	Blue/yellow	£16.10
RAC0134	HS120186	Pink/white	£16.10
RAC0136	HS120187	Blue	£16.10
RAC0138	HS120188	White	£16.10
RAC0140	HS120189	Pink	£16.10
RAC0142	HS120190	Blue	£16.10
RAC0144	HS120191	Green	£16.10

## Heathrow Scientific 4-Way Cube Tube Racks



- Cube racks, interlocking on all sides
- Rack slides together on any side allowing the user to work with different tubes on as many connected racks as desired
- Made of heavy duty polypropylene
- Autoclavable
- Dims, l x w x h: 108 x 113 x 113mm

Supplied in packs of 5

	Side 1	Side 2	Side 3	Side 4
Wells	4	10	12	16
Well Ø, mm	30	17	12.5	11.5
For Tubes, mL	50	15	3 - 7	1.5/2.0



Code	Alt Ref	Colour	Dims, w x d x h, mm	Pack	Price
RAC0469	HS29050A	Assorted	108 x 113 x 113	5	£78.00

## Heathrow Scientific 4-Way Tube Racks



- Create as large a rack as you need by linking together multiple racks
- A unique system of tabs and slots make connection easy
- Heavy duty polypropylene
- Autoclavable
- Dims, l x w x h: 174 x 95 x 52 cm

	Side 1	Side 2	Side 3	Side 4
Wells	4	12	32	32
Well Ø, mm	30	17.5	8.2	11.6
For Tubes, mL	50	15	0.5	1.5/2.0



Code	Alt Ref	Colour	Pack	Price
RAC0500	HS29022A	Natural	Each	£11.80
RAC0504	HS29022B	Blue	Each	£11.80
RAC0506	HS29022C	Green	Each	£11.80
RAC0508	HS29022F	Orange	Each	£11.80
RAC0510	HS29022D	Pink	Each	£11.80
RAC0512	HS29022E	Yellow	Each	£11.80
RAC0514	HS29022G	Assorted	5	£45.50