

SLS Select Low Form Beakers

- Ideal for general laboratory use
- White graduations and marking spot
- Made from high quality borosilicate glass
- High chemical and thermal resistance
- Staple item of any laboratory
- Standards equivalent to ISO 3819



Code	Capacity, mL	Dims, Ø x h, mm	Pack	Price
SLS5016	25	34 x 50	12	£11.61
SLS5030	50	42 x 60	12	£12.60
SLS5070	100	50 x 70	12	£14.13
SLS5072	150	60 x 80	12	£18.17
SLS5074	250	70 x 95	12	£22.30
SLS5076	400	80 x 110	12	£28.30
SLS5078	600	90 x 125	6	£17.91
SLS5080	1000	105 x 145	6	£26.50
SLS5290	5000	170 x 270	2	£158.00
SLS5292	10000	217 x 350	Each	£241.00
SLS5294	20000	285 x 430	Each	£496.00

PYREX® Low Form Beakers

PYREX®

- Manufactured from PYREX® borosilicate glass
- Complies with ISO 3819 and DIN 12331
- Excellent thermal performance and high resistance to chemical attack
- White graduations and marking spot
- Printed with Trace Code for downloadable batch certificate



Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Pack	Price
BEA12000	1000/03M	5	22 x 30	10	£43.80
BEA12020	1000/05M	10	26 x 35	10	£55.00
BEA12040	1000/11M	25	34 x 50	10	£47.00
BEA12060	1000/02M	50	42 x 60	10	£45.20
BEA12080	1000/04D	100	50 x 70	10	£44.30
BEA12100	1000/06D	150	60 x 80	10	£47.30
BEA12120	1000/10D	250	70 x 95	10	£50.00
BEA12140	1000/14D	400	80 x 110	10	£62.00
BEA12160	1000/18D	600	90 x 125	10	£76.00
BEA12180	1000/20D	800	100 x 135	10	£106.00
BEA12200	1000/22D	1000	105 x 145	10	£121.00
BEA1222	1000/28D	2000	130 x 185	Each	£21.50
BEA1224	1000/30D	3000	150 x 210	Each	£38.20
BEA1226	1000/34D	5000	170 x 270	Each	£65.00
BEA1228	1000/38D	10000	217 x 350	Each	£181.00