

Lec Laboratory Fridge-Freezer



Safe, reliable and energy efficient fridge-freezers for laboratories. These units have an anti microbial handle to help prevent the spread of bacteria and infection, come with a lock and two keys and are spark free. Other features include audible and visual high and low temperature alarms, open door alarm, reversible doors and roller castors to make the fridge-freezer easier to move for cleaning and maintenance.

- Fridge temperature range: 2 to 10°C
- Freezer temperature range: -12 to -20°C



Code	Model	Fridge Capacity, L	Freezer Capacity, L	Dims, w x d x h, mm	Price
REF4084	LSFC263UK	159	104	540 x 600 x 1680	£1077.00

Lec Pharmacy Refrigerators



Safe, reliable and energy efficient temperature controlled storage for pharmacies. Temperature can be set from 2 to 8°C set temperature - unit should operate at 6°C for optimum storage of temperature controlled drugs and vaccines.

- External digital controller
- Supplied with a lock and two keys
- Fan assisted
- Audible and visual high and low temperature alarms
- Solid door models have reversible door, glass door models do not

Essentials Range

- Solid and glass door options
- External digital controller
- Supplied with a lock and two keys
- Fan assisted
- Audible and visual high and low temperature alarms
- Reversible doors on solid door models only



Code	Model	Capacity, L	Dims, w x d x h, mm	Door Style	Price
REF4090	PEGR41UK	41	379 x 500 x 702	Glass	£1125.00
REF4001	PESR41UK	41	379 x 500 x 702	Solid	£944.00
REF3910	PESR47UK	47	470 x 465 x 530	Solid	£700.00
REF3906	PEGR82UK	82	502 x 660 x 540	Glass	£1642.00
REF3904	PESR82UK	82	502 x 660 x 540	Solid	£1173.00
REF3908	PEGR107UK	107	495 x 620 x 830	Glass	£1947.00
REF3900	PESR107UK	107	495 x 620 x 830	Solid	£1854.00