

Haier Biomedical Blood Bank Refrigerators Continued

Haier Biomedical
Intelligent Protection of Life Science

| Code | Alt Ref | Volume, L | Display | Capacity | Dims, w x d x h, mm | Shelf Positions | Price |
|---------|-----------|-----------|-------------|----------|---------------------|-----------------|-----------|
| REF2310 | HXC-149 | 149 | LED | 60 bags | 625 x 820 x 1130 | 2 | £7496.00 |
| REF2312 | HXC-149T | 149 | Touchscreen | 60 bags | 625 x 820 x 1130 | 2 | £9909.00 |
| REF2314 | HXC-429 | 429 | LED | 195 bags | 625 x 940 x 1830 | 5 | £10856.00 |
| REF2316 | HXC-429T | 429 | Touchscreen | 195 bags | 625 x 940 x 1830 | 5 | £13618.00 |
| REF2318 | HXC-629 | 629 | LED | 312 bags | 765 x 940 x 1970 | 6 | £13233.00 |
| REF2320 | HXC-629T | 629 | Touchscreen | 312 bags | 765 x 940 x 1970 | 6 | £15657.00 |
| REF2322 | HXC-1369 | 1369 | LED | 624 bags | 1545 x 940 x 1970 | 6 | £20205.00 |
| REF2324 | HXC-1369T | 1369 | Touchscreen | 624 bags | 1545 x 940 x 1970 | 6 | £23443.00 |



Haier Biomedical Pharmacy Refrigerators

Haier Biomedical
Intelligent Protection of Life Science

Small under bench pharmacy refrigerators, suitable for installation in pharmacies, pharmaceutical companies, hospitals, vaccination centres and clinics.

- Audible and visual alarms for high/low temperature, door ajar and sensor error
- Forced air cooling with optimised air distribution system designed to achieve maximum temperature uniformity and stability
- Adjustable temperature range 2 to 8°C
- Microprocessor control, LED display, solid door, and temperature adjustment with an increment of 0.1°C
- Optimised space utilisation design with multi level and adjustable shelf height to accommodate different storage requirements

| Code | Alt Ref | Capacity, L | Dims, w x d x h, mm | Price |
|---------|---------|-------------|---------------------|----------|
| SLS1130 | HYC-68 | 68 | 495 x 580 x 660 | £792.00 |
| SLS1134 | HYC-118 | 118 | 597 x 635 x 810 | £1333.00 |



Gram BIOBASIC Refrigerators

The BIOBASIC range provides cost-effective biostorage solutions for a host of applications. With forced air distribution, temperature alarms, internal and external ATEX and a user friendly interface, the BIOBASIC range offers an exceptional selection of features in a very small package.

- Temperature range: 2 to 15°C
- Acoustic door alarm, high and low temperature alarms, and a voltage free contact (NO/ NC) for remote alarms.
- Full ATEX compliance - both internal and external
- Forced air distribution system makes sure the temperature inside the cabinet remains stable and uniform
- Automatic smart defrost with re-evaporation of defrost water

| Code | Model | Capacity, L | Dims, w x d x h, mm | No. of Shelves | Price |
|---------|-------|-------------|---------------------|----------------|----------|
| REF1160 | RR210 | 125 | 595 x 640 x 831 | 2 | £1474.00 |
| REF1162 | RR310 | 218 | 595 x 640 x 1220 | 3 | £1200.00 |
| REF1164 | RR410 | 346 | 595 x 640 x 1876 | 5 | £1390.00 |