



## Liebherr Mediline Spark Free Freezers

**LIEBHERR**

Liebherr freezers with a Comfort electronic controller and spark free interior can be used for storing explosive and highly flammable substances.

- Comfort electronic controller enables temperatures to be set accurately
- ATEX has been updated. 2014 / 34 / EU (ATEX) directive, rated as II 3G Ex NA II T6
- Visual and audible alarm systems
- Integrated data memory
- Volt free alarm contact
- Static cooling
- Manual defrost
- Self closing doors with lock

Code	Model	Description	Capacity, L	Temp Range, °C	Dims, w x d x h, mm	No. of Drawers	Price
REF1720	LGUex 1500	Underbench	139	-9 to -26	601 x 618 x 820	3	£1339.00
REF1722	LGex 3410	Upright	310	-9 to -30	601 x 618 x 1840	8	£1909.00



## Liebherr Underbench Freezer

**LIEBHERR**

The Liebherr Premium underbench freezer with static cooling is an ultra-efficient unit that easily fits on or under a counter. Designed for better freezing performance and stable temperatures.

- Large surface evaporator for better cooling performance and faster freezing of the interior contents
- Highly efficient fan for consistent internal temperatures
- Adjustable temperature control, down to -26°C
- Digital temperature display with touch control buttons
- Audible alarms for high temperatures and open door
- Coated or stainless steel body with adjustable feet
- Thermoformed polystyrene liner made from a single sheet for easy cleaning
- Reversible, self closing doors have interchangeable hinges for right or left open
- Solid steel door lock for securing valuable samples and lab products

Code	Model	Capacity, L	Temp Range, °C	Dims, w x d x h, mm	Price
REF1730	GGU 1500 Premium	143	-9 to -26	600 x 615 x 830	£879.00



## Liebherr Laboratory Upright Freezers with Static Cooling

**LIEBHERR**

Liebherr freezers offer perfect cooling performance and meet the most exacting demands. The interior temperature can be precisely set to between -14 and -28°C. Premium quality materials, optimised and high performance components, and first rate workmanship ensure a long service life and ensure that appliances offer optimum storage conditions for frozen products.

- Integrated temperature alarm with both acoustic and visual signals
- Static cooling
- Manual defrost
- Factory fitted with door hinges on the right, but with the ability to change the hinge side
- SwingLine design with elegantly sweeping lines, ergonomically designed, easy to clean, and easy to open with full-length integrated grip handle
- Electronic control system with external digital temperature display
- Environmentally friendly R600a or R290 refrigerants
- Capacity: 513L