



Gram BIOBASIC Freezers

The BIOBASIC range provides cost effective biostorage solutions for a host of applications. With forced air distribution, temperature alarms, internal and external ATEX and a user friendly interface, the BIOBASIC range offers an exceptional selection of features in a very small package.

- Acoustic door alarm, high and low temperature alarms, and a voltage free contact (NO/NC) for remote alarms
- Full ATEX compliance – both internal and external
- Forced air distribution system makes sure the temperature inside the cabinet remains stable and uniform
- Automatic smart defrost with re-evaporation of defrost water
- Temperature range: -5 to -25°C

Code	Model	Capacity, L	Dims, w x d x h, mm	No. of Shelves	Price
REF1166	RF210	125	595 x 640 x 831	2	£1540.00
REF1168	RF310	218	595 x 640 x 1220	3	£1680.00
REF1170	RF410	346	595 x 640 x 1876	5	£1820.00



Labcold Basic Underbench Laboratory Freezer



Built for value and ideal for when there is no risk that volatile, flammable materials or components might be stored. Engineered to provide good temperature control, this model provides basic frozen storage at a lower price while still providing the confidence of the Labcold name.

- Enamelled steel exterior and plastic interior
- Robust construction
- Operator controlled manual defrost
- 2 clear plastic drawers and top compartment with drop down plastic flap
- Two adjustable feet at the front

Code	Capacity, L	Temp Range, °C	No. of Drawers	Dims, w x d x h, mm	Price
RLVL03203	70	-18 to -23	2 + 1 flap	500 x 490 x 835	£511.00



Labcold Basic Upright Laboratory Freezer



Built for value and ideal for when there is no risk that volatile, flammable materials or components might be stored. Engineered to provide good temperature control, this model provides basic frozen storage at a lower price without significant sacrifices in the important areas of function and reliability.

- Robust steel exterior, easy clean plastic interior
- Operator controlled manual defrost
- Six fixed shelves
- Door lock with two keys
- Digital temperature display
- High/low temperature alarm, door open alarm
- 4 swivel wheels, 2 with brakes
- Solid, self closing door

Code	Capacity, L	Temp Range, °C	No. of Shelves	Dims, w x d x h, mm	Price
RLVF2025	543	-18 to -24	6	750 x 755 x 1960	£2282.00