

## Haier -40°C Biomedical Freezers

Haier Biomedical  
Intelligent Protection of Life Science

Haier Biomedical UK's range of -40°C biomedical freezers offer a large capacity storage space with excellent temperature uniformity. The freezers utilise shelf integrated evaporators to offer rapid pull down, recovery time and ensure your samples are stored at the correct temperature no matter where they are within the freezer. This range is designed to store vaccines, blood plasma and many other biological materials.

- Temperature range: -20 to -40°C
- Rapid cooling with integrated evaporator shelves
- Improved temperature uniformity with dual capillary tube design
- Improved door seal design with 2 gaskets reduces frost build
- LED digital display
- Optional USB data logging on some models
- Multiple function alarms include high/low temp, sensor error (all models)
- Larger models also have low battery, high ambient temperature, power failure alarms



Code	Alt Ref	Description	Capacity, L	USB Datalogging	Dims, w x d x h, mm	Price
SLS1122	DW-40L92	Underbench	92	No	640 x 610 x 860	£1482.00
SLS1124	DW-40L262	Upright	262	No	700 x 705 x 1665	£2358.00
SLS1126	DW-40L348	Upright double door	348	Yes	765 x 810 x 1860	£5256.00

## Eppendorf CryoCube® 570/660 Series

eppendorf

To protect your precious samples, CryoCube® freezers are built with reliability in mind. Constructed from premium components, they ensure consistent performance and preservation year after year. The CryoCube® F570h/hw and FC660h (high energy efficiency models) further improve energy efficiency by up to 25%, compared to the F570n and FC660. They combine all the enhanced technology of the F570n and the FC660 with additional vacuum insulation coverage.

- Quick door opening and improved temperature uniformity with new automatic vent port located within the user interface
- Reduced energy consumption with high efficiency fan, compressor, and condenser
- Ergonomic door handle of CryoCube® models based on Eppendorf PhysioCare Concept® requires less force and provides easier access
- Easy access to samples with new magnetic closures on inner doors
- Green cooling (R290/R170) and green foam insulation for high levels of sustainability
- Flexible shelf positions for customisable setup



**Models featuring optional accessory features such as chart recorders, CO<sub>2</sub> back up systems and liquid nitrogen back-up systems are available – see SLS website for full range listing.**

Code	Alt Ref	Model	Style	Capacity, L	Compartments	Max Racks per Freezer	Energy Consumption, kWh/24hr	Dims, w x d x h, mm	Price
FRE6200	F571300031	F570n	Upright	570	5	25	8.3	1103 x 898 x 1950	£12524.00
FRE6214	F571340031	F570h	Upright, high efficiency	570	5	25	7.4	1103 x 898 x 1950	£13579.00
FRE6216	F571340131	F570hw	Upright, high efficiency, water cooled	570	5	25	7.8	1103 x 898 x 1950	£14696.00
FRE6058	F660300001	FC660	Chest	660	—	40	10.6	2070 x 845 x 1090	£13336.00
FRE6060	F660320001	FC660h	Chest, high efficiency	660	—	40	9.0	2070 x 845 x 1090	£14918.00