



A

DURAN® Low Form Beakers



- Easy filling is possible via the large vessel opening
- Complies with ISO 3819
- Manufactured from borosilicate glass 3.3
- Spout for clean pouring
- With easy to read graduations and large labelling field for easy marking in fired on, highly durable, white ceramic
- Retrace code which allows the downloading of an individual batch certificate

Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Pack	Price
BEA11000	211060701	5*	22 x 30	10	£50.50
BEA11020	211060804	10*	26 x 35	10	£58.00
BEA11040	211061406	25	34 x 50	10	£49.00
BEA11060	211061706	50	42 x 60	10	£46.90
BEA11080	211062402	100	50 x 70	10	£46.60
BEA11100	211062908	150	60 x 80	10	£49.00
BEA11120	211063604	250	70 x 95	10	£48.70
BEA11140	211064103	400	80 x 110	10	£63.00
BEA11160	211064806	600	90 x 125	10	£77.00
BEA11180	211065305	800	100 x 135	10	£105.00
BEA11200	211065408	1000	105 x 145	10	£121.00
BEA1122	211066301	2000	132 x 185	Each	£21.00
BEA1124	211066807	3000	152 x 210	Each	£37.90
BEA1126	211067306	5000	170 x 270	Each	£67.00
BEA1128	211068602	10000	217 x 350	Each	£193.00

* Without graduation or retrace code



DURAN® Silicone Lids



- For covering and sealing laboratory vessels
- Reusability makes them an environmentally friendly, cost effective alternative to one time use covers and films
- Stretchable and thus suitable for a large number of vessel openings
- Slightly larger tab can be used as a writing field
- Available in different colours for clear differentiation and identification
- Easy to remove with the aid of a flap, both with and without gloves
- Resistant to both chemicals and heat
- Recommended for use at temperatures between -40 and +180°C
- Autoclavable (silicone lid only without glassware)
- Manufactured from silicone

Code	Alt Ref	Size	Colour	Dims, Ø x h, mm	Price
BOT1248	291101119	Small	Pink	22.5 x 42.2	£1.70
BOT1258	291101127	Small	Cyan	22.5 x 42.2	£1.70
BOT1260	291101135	Small	Green	22.5 x 42.2	£1.70
BOT1262	291102115	Medium	Pink	63.7 x 22.5	£2.20
BOT1264	291102123	Medium	Cyan	63.7 x 22.5	£2.20
BOT1266	291102131	Medium	Green	63.7 x 22.5	£2.20
BOT1268	291103111	Large	Pink	83.7 x 23.2	£2.70
BOT1270	291103128	Large	Cyan	83.7 x 23.2	£2.70
BOT1272	291103136	Large	Green	83.7 x 23.2	£2.70