



SLS Select Nitrile Gloves



- Blue nitrile
- Powder and latex free
- Good resistance to chemical substances
- Protection from infections
- Suitable for applications involving certain cytostatic agents
- 4.8g in weight
- AQL 1.5
- EN 455 (medical gloves for single use)
- EN 374

Code	Size	Pack	Price
SAF8578	Small	100	£20.28
SAF8580	Medium	100	£20.28
SAF8582	Large	100	£20.28
SAF8584	X large	90	£20.28



SLS Lab Pro Nitrile Gloves



- Blue nitrile
- Powder and latex free
- Good resistance to chemical substances
- Protection from infections
- Suitable for applications of certain cytostatic agents
- 4.8g in weight
- AQL 0.65
- EN 455 (medical gloves for single use)
- EN 374

Code	Size	Pack	Price
SAF8586	Small	100	£22.00
SAF8588	Medium	100	£22.00
SAF8590	Large	100	£22.00
SAF8592	X large	90	£22.00



Ansell TouchNTuff® 92-600 Nitrile Gloves



- Teal green nitrile, double chlorinated
- TNT™ technology - a proprietary polymer formulation that provides superior splash resistance for durable protection with a soft, comfortable feel
- Silicone, latex and powder free
- Beaded cuff
- Ambidextrous
- Superior comfort and durability
- Anti static
- Length: 240mm
- Category III
- AQL 1.5
- EN ISO 374-1:2016 Type B JKPT (chemical risks)
- EN ISO 374-5:2016 VIRUS (protection from viruses, bacteria and fungi)
- Food contact approved

Supplied in packs of 100

Code	Alt Ref	Size	Price
SAF1114	92-600-S	Small (6.5 - 7)	£66.00
SAF1116	92-600-M	Medium (7.5 - 8)	£66.00
SAF1118	92-600-L	Large (8.5 - 9)	£66.00
SAF1120	92-600-XL	X large (9.5 - 10)	£66.00