

## Glove Plus® Ultra Nitrile Gloves



- High performance protection against a wide range of chemicals and blood borne pathogens
- Durable and strong
- Textured fingertips for enhanced grip
- Below 2mg/g glove powder content
- Enhanced processing, silicon free, ultra performance
- Suitable for contact with chemicals and blood borne pathogens
- Medium/high risk category
- Duration of use: 120+ minutes
- Length: 240mm
- PPE class III
- EN 455 parts 1 & 2
- AQL 0.65
- ISO 16604b:2003 viral test passed
- EN 420 dexterity test passed level 5
- ASTM D6978-05 permeation of chemotherapy drugs passed



Supplied in packs of 100

| Code    | Alt Ref | Size    | Price  |
|---------|---------|---------|--------|
| SAF6322 | UNP32   | Small   | £41.00 |
| SAF6324 | UNP33   | Medium  | £41.00 |
| SAF6326 | UNP34   | Large   | £41.00 |
| SAF6328 | UNP35   | X large | £41.00 |

## Glove Plus® Lite Nitrile Gloves



- The perfect medium weight nitrile gloves
- Superb chemical barrier properties and precision grip, they are the ideal everyday glove.
- Below 2mg/g glove powder content
- Textured fingertips and precision grip
- Medium risk category
- Suitable for contact with chemicals and blood borne pathogens
- Duration of use: 30 – 120 minutes
- Length: 240mm
- Medical grade
- EN 455 parts 1 & 2
- EU PPE class III
- AQL 0.65
- ISO 16604b:2003 viral test passed
- EN 420 dexterity test passed level 5
- ASTM D6978-05 permeation of chemotherapy drugs passed



Supplied in packs of 100

| Code    | Alt Ref | Size    | Price  |
|---------|---------|---------|--------|
| SAF5682 | LN302   | Small   | £38.90 |
| SAF5684 | LN303   | Medium  | £38.90 |
| SAF5686 | LN304   | Large   | £38.90 |
| SAF5688 | LN305   | X large | £38.90 |

## Glove Plus® Prime Nitrile Gloves



- The ideal glove for use over short periods of time where splash contact from aggressive chemicals is unlikely
- Barrier protection against blood borne pathogens
- Below 2mg/g glove powder content
- Increased sensitivity, fingertip textured for precision grip
- Duration of use: 0 – 30 minutes
- Suitable for contact with chemicals and blood borne pathogens
- Low risk category
- Medical grade
- Length: 240mm
- EN 455 parts 1 & 2
- EU PPE class III
- AQL 0.65
- ISO 16604b:2003 viral test passed
- EN 420 dexterity test passed level 5



Supplied in packs of 100

| Code    | Alt Ref | Size    | Price  |
|---------|---------|---------|--------|
| SAF5520 | PN301   | X small | £18.90 |
| SAF5522 | PN302   | Small   | £18.90 |
| SAF5524 | PN303   | Medium  | £18.90 |
| SAF5526 | PN304   | Large   | £18.90 |
| SAF5529 | PN305   | X large | £18.90 |