



Ebony Black Nitrile Gloves

- High performance, high quality matt black nitrile glove
- Great for a range of purposes due to its versatility, durability and sensitivity
- Gully textured finish provides enhanced grip
- 100% NR latex free
- Below 2mg/g glove powder content
- Fully textured, premium content
- Duration of use: 120+ minutes
- Suitable for contact with chemicals and blood borne pathogens
- Medium/high risk category
- Length: 240mm
- EN 455 parts 1 & 2
- AQL 0.65
- EU PPE class III
- ISO 16604b:2003 viral test passed
- EN 420 dexterity test passed level 5

Supplied in packs of 100

Code	Alt Ref	Description	Price
SAF5404	EB702	Small	£41.20
SAF5406	EB703	Medium	£41.20
SAF5407	EB704	Large	£41.20
SAF5408	EB705	X large	£41.20



SHIELDskin™ Chem™ NEO NITRILE™ 300 Gloves



- Red polychloroprene/nitrile
- Powder and latex free
- Textured fingertips for grip
- Beaded cuff
- Ambidextrous
- twinSHIELD™ double walled protection
- Length: 300mm
- Category III
- AQL 0.25
- EN ISO 374-1:2016 + A1:2018 Type A AJKLNPT (chemical risks)
- EN 16523-1:2015 + A1:2018 (chemical permeation)
- ISO 374-4:2019 (resistance to degradation by chemicals)
- EN ISO 374-5:2016 VIRUS (protection from viruses, bacteria and fungi)
- ISO 16604:2004 Procedure B (viral penetration test)
- EN 388:2016 2000X (mechanical risks)
- EN 455-1/2/3/4:2000 (medical gloves for single use)

Supplied in packs of 40

Code	Alt Ref	Size	Price
SAF1068	669251	X small [6]	£110.00
SAF1070	669252	Small [7]	£110.00
SAF1072	669253	Medium [8]	£110.00
SAF1074	669254	Large [9]	£110.00
SAF1076	669255	X large [10]	£110.00
SAF1078	669256	2X large [11]	£110.00