


## duoSHIELD™ PFT Nitrile 240 Gloves



- Blue nitrile
- Powder and latex free
- Textured fingertips for grip
- Beaded cuff
- Ambidextrous
- Length: 240mm
- Category III
- AQL 0.65
- EN ISO 374-1:2016 + A1:2018 Type B KPT (chemical risks)
- EN 16523-1:2015 + A1:2018 (chemical permeation)
- ISO 374-4:2019 (resistance to degradation by chemicals)
- EN ISO 374-5:2016 VIRUS (protection from viruses, bacteria and fungi)
- ISO 16604:2004 Procedure B (viral penetration test)
- EN 455-1/2/3/4:2000 (medical gloves for single use)



Supplied in packs of 200

| Code           | Alt Ref | Size         | Price  |
|----------------|---------|--------------|--------|
| <b>SAF3102</b> | 658121  | X small (6)  | £75.00 |
| <b>SAF3104</b> | 658122  | Small (7)    | £75.00 |
| <b>SAF3106</b> | 658123  | Medium (8)   | £75.00 |
| <b>SAF3108</b> | 658124  | Large (9)    | £75.00 |
| <b>SAF3114</b> | 658125  | X large (10) | £75.00 |

## Kimtech™ Purple Nitrile™ Gloves

**KIMTECH™**

- Purple synthetic nitrile polymer
- Powder, latex and silicone free
- Textured fingertips for enhanced grip
- Beaded cuff
- Ambidextrous
- Anti static
- Length: 240mm (XS - M), 250mm (L - XL)
- Category III
- AQL 0.65
- EN ISO 374-1:2016 Type B JKT (chemical risks)
- EN 374-4:2014 (resistance to degradation by chemicals)
- EN ISO 374-5:2016 VIRUS (protection from viruses, bacteria and fungi)
- Food contact approved



| Code           | Alt Ref | Size    | Pack | Price  |
|----------------|---------|---------|------|--------|
| <b>SAF3408</b> | 90625   | X small | 100  | £64.00 |
| <b>SAF3410</b> | 90626   | Small   | 100  | £64.00 |
| <b>SAF3412</b> | 90627   | Medium  | 100  | £64.00 |
| <b>SAF3414</b> | 90628   | Large   | 100  | £64.00 |
| <b>SAF3416</b> | 90629   | X large | 90   | £64.00 |

## Electrophoresis

A complete range of electrophoresis tanks, power supplies and accessories brought to you by SLS.

See pages **341-382**

