

Dermagrip Nitrile Gloves

- An alternative solution for individuals who are allergic to natural rubber latex
- Formulated with 100% nitrile (acrylonitrile-butadiene)
- Do not contain natural latex protein
- Do not exhibit any chemical residue level when tested with HPLC using water extraction method
- Biocompatibility friendly in addition to providing the usual barrier protection property needed for examination gloves
- Ambidextrous
- Textured finger surface
- Beaded cuffs

Supplied in packs of 100




Code	Alt Ref	Description	Price
SAF9366	GL91250S	Small	£26.00
SAF9368	GL91250M	Medium	£26.00
SAF9370	GL91250L	Large	£26.00

Handsafe GN830 Long Cuff Nitrile Gloves

- Manufactured from nitrile butadiene
- Accelerator free, removing the risk of contact dermatitis caused by these chemicals
- Latex free, reducing the risk of hypersensitivity and allergic reactions
- Suitable for use with cytotoxic/chemotherapy drugs
- Powder free to reduce the risk of dust contamination
- Textured finish for improved grip
- Long cuff: 29cm for extra protection
- PPE category III
- Compliant with the Medical Devices Regulation (EU) 2017/745 class 1
- AQL 1.5
- EN 455: Medical gloves for single use parts 1, 2, 3 & 4
- EN 420:2003+A1:2009
- EN ISO 374-1:2016 Type C K (chemical risks)
- EN ISO 374-5:2016 VIRUS (protection from viruses, bacteria and fungi)
- EN 1186 - suitable for food contact
- Product food approved to EC1935/2004
- EU type examination carried out by SATRA

Supplied in packs of 100



Code	Size	Price
SAF8608	Small	£32.00
SAF8610	Medium	£32.00
SAF8612	Large	£32.00
SAF8614	X large	£32.00

KEY

Throughout the catalogue the following symbols are used:



New product



Service and repairs available



Calibration service available



Made in Britain



Sustainable features



Autoclavable



No. years warranty



No. years warranty can be extended



Lifetime warranty