




A

SLS Select Beakers



- Polypropylene
- Autoclavable

Code	Capacity, mL	Pack	Price
BEA1060	50	12	£7.40
BEA1062	100	12	£8.80
BEA1064	250	12	£12.90
BEA1067	500	4	£8.41
BEA1068	1000	4	£15.80
BEA1070	2000	4	£28.70
BEA1072	5000	2	£38.80
BEA1074	10000	Each	£38.60



A

Azlon® Beakers, Tapered



- Tapered shape to make efficient use of your storage space
- Comply to ISO 7056
- Autoclavable at 121°C with good chemical resistance
- Stackable space saving design
- Made from high clarity polypropylene
- Durable and economical with clear, moulded permanent graduations
- Precise non drip pouring spouts

Code	Alt Ref	Capacity, mL	Graduations, mL	Dims, Ø x h, mm	Pack	Price
BEA5120	BPM0025P	25	5	41.5 x 51	10	£9.30
BEA5122	BPM0050P	50	10	50.5 x 61	10	£11.50
BEA5124	BPM0100P	100	20	61.5 x 73	10	£14.60
BEA5126	BPM0250P	250	50	80.5 x 96	10	£21.10
BEA5128	BPM0400P	400	50	92.5 x 111	10	£29.20
BEA5130	BPM0500P	500	50	102.5 x 121	10	£33.10
BEA5132	BPM0600P	600	50	104 x 126	5	£19.00
BEA5134	BPM1000P	1000	100	120 x 148	5	£21.90
BEA5136	BPM2000P	2000	200	163 x 188	Each	£8.10

Authorised distributor



We offer a broad portfolio of products from key suppliers such as:



Honeywell

MERCK

Sigma-Aldrich
Lab & Production Materials

Supelco
Analytical Products

spex
an antylla scientific company