

Dupont™ Tyvek® 500 Lab Coats



- Flashspun high density polyethylene
- Permeable to both air and water vapour, yet repels water based liquids and aerosols
- Excellent barrier against fine particles and fibres
- Ultra-low linting
- Anti statically treated
- Press stud closures
- 3 pockets
- Collar
- Stitched internal seams
- Category III
- EN 14126:2003 (infective agents)
- EN 13034:2005 + A1:2009 Type 6 (PB) (limited chemical splash protection)
- EN 1149-1 (antistatic treatment)



Supplied in packs of 50

Code	Alt Ref	Size	Price
SAF0302	D13396079	Medium	£416.00
SAF0304	D13396069	Large	£416.00
SAF0306	D13396050	X large	£416.00

Unisex Howie Lab Coats

- Polyester/cotton mix
- Weight: 245gsm
- One jetted chest breast pocket and two jetted hip pockets
- Offset nine stud fastening with broad underlay
- Knitted cuffs
- Back darts
- Other colours available – please enquire



Regular

- Length: 1067mm

Code	Size	Chest Size, cm	Pack	Price
SAF0037	X small	76	Each	£49.00
SAF0038	X small	84	Each	£49.00
SAF0040	Small	88	Each	£49.00
SAF0042	Small	92	Each	£49.00
SAF0044	Medium	96	Each	£49.00
SAF0046	Medium	100	Each	£49.00
SAF0048	Large	108	Each	£49.00
SAF0050	Large	112	Each	£49.00
SAF0052	X large	116	Each	£49.00
SAF0054	X large	124	Each	£49.00
SAF0055	2X large	128	Each	£49.00
SAF0056	2X large	132	Each	£49.00
SAF0058	3X large	140	Each	£49.00

Tall

- Length: 1118mm

Code	Size	Chest Size, cm	Pack	Price
SAF2660	X small	76	Each	£25.00
SAF2662	X small	84	Each	£25.00
SAF2664	Small	88	Each	£25.00
SAF2666	Small	92	Each	£25.00
SAF2668	Medium	96	Each	£25.00
SAF2670	Medium	100	Each	£25.00
SAF2672	Large	108	Each	£25.00
SAF2674	Large	112	Each	£25.00
SAF2676	X large	116	Each	£25.00
SAF2678	X large	124	Each	£25.00
SAF2680	2X large	128	Each	£25.00
SAF2682	2X large	132	Each	£25.00
SAF2684	3X large	140	Each	£25.00