

## Asecos SLX Corrosive and Alkali Recirculating Filtration Cabinet

*asecos*

The SLX range offers superior storage and filtration of corrosive liquids and alkalis. The liquids should be segregated in accordance with HSG71 and the safety data sheet.

- Choice of 3 widths: 1197 and 600mm, all depths 615mm and all heights 2300mm
- Two separate hermetically sealed cabinets in one
- Low metallic content
- Pull out tray shelves with melamine drip tray
- Light grey cabinet
- Multi stage wideband filter, touch display monitor for service and user operation settings, air and fan monitoring, filter change required annually
- Ideal for where no on site extraction is available
- Delivery service – to goods in or point of use where the cabinet is delivered to where it is required, it is levelled and the packaging removed (site survey required for point of use)



Please enquire for more options

Code	Alt Ref	Description	Dims, w x d x h, mm	Standards	Price
SAF6038	32939-001-30615	SLX recirculating corrosive and alkali cabinet 600mm widthh with 3 pull out tray shelves in the top and 3 in the bottom complete with filtration unit	600 x 615 x 2300	CE	£6274.00
SAF6040	32953	Replacement filter - activated carbon	—	—	£798.00

## Asecos Under Bench 90 Minute Fire Resistant Cabinet

*asecos*

The UB range offers a broad range of under bench storage for flammable liquids, it is manufactured to BS EN 14470 part 1, the cabinets are ideal for storage where space is premium. They have proven themselves in many modern laboratory designs, whether building a new laboratory or refurbishing an existing laboratory they are the perfect solution. Contents must be stored in line with HSG71 and the chemical safety data sheets.

- 90 minute fire resistance
- Choice of widths: 1402, 1102, 893 and 593mm, choice of depths 500 and 574mm, standard height 600mm, 800mm height available only in 1102 and 593mm widths, also with wing doors or pull out draw front
- Ready for connection to an extraction system through a 50mm port
- Self closing doors, self damping and self sealing in the event of fire
- Light grey cabinets
- Internal options of second level draws and fixed shelves
- All cabinets contain a sump and perforated insert avoiding displacement of sump volume in the event of a spill
- Delivery service – to goods in or point of use where the cabinet is delivered to where it is required, it is levelled and the packaging removed (site survey required for point of use)



Please enquire for more options

Code	Alt Ref	Description	Dims, w x d x h, mm	Standards	Price
SAF6042	32658-001-35378	Under bench cabinet with pull out draw front and second level drawer with perforated insert	593 x 600 x 800	BS EN 14470 part 1	£2437.00