



Asecos UFA Recirculating Extraction for Flammable Liquids

asecos

The UFA is a recirculating extraction for tall flammable storage and under bench cabinets, a compact solution and cost effective should you not have access to on-site extraction.

- Compact design
- Choice of two options, one for tall cabinets and one for under bench cabinets
- Activated carbon filtration for flammable atmospheres only
- Quiet in operation only 39dBA
- Constant filter and fan monitoring with both visual and audible warning
- Light grey coloured unit
- Delivery service – to goods in

Code	Alt Ref	Description	Dims, w x d x h, mm	ATEX II Compliant	Price
SAF6044	25099	UFA 25099 recirculating air filter system with exhaust air monitoring, for tall cabinets	305 x 555 x 210	Yes	£1854.00
SAF6046	27144	UFA 27144 recirculating air filter system with exhaust air monitoring, for under bench cabinets	305 x 555 x 210	Yes	£1895.00
SAF6048	24121	UFA replacement filter (activated carbon) with reset cable and manual	—	—	£204.00



Asecos CS Classic Cabinet for Non-Flammable and Non-Corrosive Toxic Chemicals

asecos

The CS Classic range is designed to offer a wide range for the storage of chemicals namely toxics that are non flammable or non corrosive. The range offers ventilation and can be connected to extraction and also offer load tested shelving.

- Choice of 3 widths: 1055, 810 and 545mm, available heights 1105 and 1950mm, all cabinets in this range are 520mm deep
- Choice of 7 door colours, anthracite grey surround
- Can be easily relocated within a building
- Packages with fixed shelves or drawers are available in powder coated sheet steel
- Options include bottle racks, drawer mats and internal storage boxes
- Delivery service – to goods in or point of use where the cabinet is delivered to where it is required, it is levelled and the packaging removed (site survey required for point of use)

Please enquire for more options

Code	Alt Ref	Description	Dims, w x d x h, mm	Standards	Price
SAF6062	30495-041-30508	CS Classic cabinet with 3 fixed shelves, bottom collecting sump and perforated insert	1055 x 520 x 1950	CE	£1415.00

KEY

Throughout the catalogue the following symbols are used:



New product



Service and repairs available



Calibration service available



Made in Britain



Sustainable features



Autoclavable



No. years warranty



No. years warranty can be extended



Lifetime warranty