



A

Thermo Scientific™ Nalgene™ Low Form Beakers, Griffin Style

Thermo
SCIENTIFIC

- Single spout design
- Silk screened with graduations, size code, resin code, catalog number, max. temperature and "No Flame" symbol for easy identification and safe use
- Rugged, reusable, construction
- High quality, laboratory grade plastic materials for dependably low leachables and extractables
- Meets ISO/DIN 7056 International Standards for Laboratory Plasticware Accuracy
- Ribs for easy stacking
- Autoclavable

Code	Alt Ref	Capacity, mL	Graduations, mL	Pack	Price
BEA3340	1201-0030	30	5	12	£36.20
BEA3342	1201-0050	50	5	12	£39.10
BEA3344	1201-0100	100	10	12	£38.20
BEA3346	1201-0150	150	10	12	£46.90
BEA3348	1201-0250	250	25	6	£30.50
BEA3350	1201-0400	400	25	6	£38.90
BEA3352	1201-0600	600	50	4	£33.90
BEA3354	1201-1000	1000	50	3	£32.00
BEA3356	1201-2000	2000	100	Each	£23.30
BEA3358	1201-4000	4000	250	Each	£40.40



A

Simport® Scientific Tricorn™ Beakers

Simport
Scientific

- Provide three dripless pouring spouts
- Tough, and unbreakable
- Suitable for use with commonly used acids, alkalis and solvents
- Autoclavable
- Made of translucent polypropylene

Supplied in packs of 100

Code	Alt Ref	Capacity, mL	Graduations, mL	Dims, Ø x h, mm	Price
BEA5000	B700-50	50	5	49 x 57	£50.00
BEA5002	B700-100	100	10	58 x 72	£71.00
BEA5004	B700-250	250	10	76 x 96	£96.00
BEA5006	B700-400	400	20	85 x 108	£96.00
BEA5008	B700-800	800	50	107 x 133	£134.00
BEA5010	B700-1L	1000	50	115 x 145	£143.00



A

Nickel Electro Metal Beakers

NICKEL-ELECTRO LTD.

- Stainless steel 304 spec 18/8 non magnetic, bright finish
- Durable seamless design with a flat base for stability
- Even wall thickness and generous flared rim for durability
- Spout for easy pouring
- Will not shatter or chip
- All surfaces are non reactive and offer zero absorption
- Can be sterilised by autoclaving

Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Price
BEA3320	201	200	63 x 73	£45.00
BEA3322	202	350	65 x 115	£48.90
BEA3324	203	550	77 x 124	£64.00
BEA3326	204	700	75 x 172	£58.00
BEA3328	206	1000	110 x 135	£92.00
BEA3330	205	1250	100 x 162	£94.00
BEA3332	207	2000	126 x 156	£117.00
BEA3336	209	5000	177 x 206	£163.00