

SHARPSGUARD® Sharps Containers

- Enhanced to enable compliance with the requirements of the Safe Management of Healthcare Waste guidance published by the Department of Health
- Colour coded containers, lids and labels are available in a variety of sizes, whilst retaining the very highest standards
- Manufactured to BS EN ISO 23907 and UN 3291

Yellow

- For the disposal of sharps that have been contaminated with medicines (not cytotoxic or cytostatic)
- For disposal by incineration or other authorised process

Code	Alt Ref	Capacity, L	Dims, h x Ø, mm	Pack	Price
SAF7272	DD472YL	2.5	166 x 176	48	£102.00
SAF7273	DD471YL	5	285 x 176	48	£147.00
SAF7274	DD473YL	7	207 x 248	40	£145.00
SAF7275	DD476YL	11.5	249 x 275	20	£102.00
SAF7277	DD475YL	22	343 x 327	10	£100.00



Orange

- For the disposal of sharps that have not been contaminated with medicines
- For disposal by incineration or other authorised process
- Com-plus is rectangular for easy storage in medical bag and has a sliding door with an aperture which will take a 20mL syringe
- Theatre+ is designed to be used upright to ease the disposal of long instruments such as staple guns as well as chest drains, it is ideal for theatres, acute medical areas and x-ray departments

Code	Alt Ref	Capacity, L	Dims, mm	Pack	Price
SAF7352	DD477OL	1	w x d x h: 68 x 204 x 119	Each	£2.40
DD4790	DD479OL	1 com-plus	w x d x h: 68 x 204 x 120	30	£47.30
SAF72720	DD472OL	2.5	h x Ø: 166 x 176	48	£103.00
DD4740L	DD474OL	3.75	h x Ø: 221 x 176	48	£95.00
SAF72730	DD471OL	5	h x Ø: 285 x 176	48	£147.00
DD4730L	DD473OL	7	h x Ø: 207 x 248	40	£116.00
DD4760L	DD476OL	11.5	h x Ø: 249 x 275	20	£81.00
DD4750L	DD475OL	22	h x Ø: 343 x 327	10	£83.00
SAF7004	DD5400L	32 theatre+	w x d x h: 291 x 307 x 710	4	£106.00



Sharpak® Sharps Containers, Orange

- High quality regulatory sharps containers for single use, safe collection and disposal
- For fully discharged sharps not contaminated with cytotoxic and cytostatic medicines or sharps not contaminated with medicinal products
- For disposal of by incineration or other approved method
- VisiLOCK simple lock indicator system
- One piece base and handle moulding
- Quick and easy to assemble
- Sequential numbering for traceability
- Approved to ISO 23907:2012
- Tested and certified in accordance with the provisions of the United Nations recommendations in the Transportation of Dangerous Goods, Chapter 9, UN3291 Clinical Waste, Unspecified, NOS

Code	Alt Ref	Capacity, L	Capacity to Fill Line, L	Dims, w x d x h, mm	Price
SAF7114	SHK/CA1526SCO	6	4.5	220 x 234 x 229	£7.09
SAF7116	SHK/CA1536SCO	12	9.6	258 x 325 x 268	£9.56
SAF7118	SHK/CA1548SCO	22	17.6	315 x 327 x 383	£18.48

