



Sampling Systems SteriWare® Volumetric Spoons



- Enable set volumes to be accurately measured and sampled
- Suitable for use with powders, granules, liquids, gels and suspensions
- Manufactured from high quality FDA, EC and food grade polymers
- Long handle keeps hands away from the material being sampled or measured
- Designed to be used straight from the bag which eliminates the risk of cross contamination and means there is no costly cleaning
- Manufactured and packed in a cleanroom
- Individually packaged

Supplied in packs of 100

Code	Alt Ref	Volume, mL	Sterility	Price
SAM1082	8640H-01	1	Non sterile	£53.00
SAM1083	8640H-01S	1	Sterile	£66.00
SAM1084	8640H-0123S	1.23	Sterile	£66.00
SAM1085	8640H-0246S	2.46	Sterile	£66.00
SAM1086	8640H-05	5	Non sterile	£53.00
SAM1087	8640H-05S	5	Sterile	£66.00
SAM1088	8640H-10	10	Non sterile	£53.00
SAM1089	8640H-10S	10	Sterile	£66.00
SAM1090	8640H-20	20	Non sterile	£53.00
SAM1091	8640H-20S	20	Sterile	£66.00



Azlon® Polypropylene Scoops



- White polypropylene
- Flat base for horizontal standing and a flat top for levelling off powders
- Capacity is clearly moulded on each size

Code	Capacity, mL	Overall Length, mm	Pack	Price
SCO1004	10	102	5	£6.00
SCO1006	25	137	5	£10.00
SCO1008	50	160	5	£12.00
SCO1010	100	200	Each	£5.00
SCO1012	250	260	Each	£8.00
SCO1014	500	315	Each	£13.00
SCO1016	1000	385	Each	£22.00



Stainless Steel Scoops

- Light duty 304 grade stainless steel
- Cylindrical

Code	Capacity, mL	Price
SCO2000	100	£30.00
SCO2002	250	£35.00
SCO2004	500	£38.00
SCO2006	1000	£47.00